

PN3070 MALS

Multi-Angle Light Scattering Detector



PN3070 MALS Detector

Specifications

- Sample Cell
100 μ L internal volume, PEEK cell with max. pressure up to 35bar / 500psi
- Static Light Scattering Volume
20 nL
- Detection Principle
Multi-Angle Static Laser Light Scattering Detection
- Light Scattering Angles
35°, 50°, 75°, 90°, 105°, 130°, 145°
- Precision of Molar Mass
+/- 1 %
- Molar Mass Range
10³ to 10⁹ Daltons
- Radius of Gyration Range
10 nm to 500 nm
- Hydrodynamic Radius Range
Rh: 13 nm to 195 nm
Dh: 26 nm to 390 nm
- Laser Specifications
35 mW, 635 nm
- Laser Life Time
Approx. 9.000 hours
- Solvent Compatibility
All commonly used aqueous/organic solvents
- Electronics
4 standard inputs with 16-24 bit resolution and 16 optional analog channels with 24 bit resolution. Software selectable gain adjustment! Suitable for most commonly used detectors as Refractive Index, UV/VIS and Viscometer, etc..
- Total Size
ca. 21 cm x 44 cm x 11 cm
- Shipping Weight
ca. 8,5 kg
- Power requirement
230/110 V, 50/60 Hz
- Available Options
- Temperature Control 50°C +/-0.015°C
- Expansion Board: 15 analog inputs with 24 bit
- Requires a PC with USB port for operation!



Postnova Analytics GmbH
Max-Planck-Str. 14
86899 Landsberg, GER
T: +49 8191 985 688 0
F: +49 8191 985 688 99

Postnova Analytics Inc.
230 S, 500 E, Suite 120
84102 Salt Lake City, USA
T: +1 801 521 2004
F: +1 801 521 2884


info@postnova.com
www.postnova.com

The PN3070 MALS Static Laser Light Scattering Detector is a new type of detection system, which combines high sensitivity fiber optics with advanced and unique CCD detection. The system enables the simultaneous and independent measurement of up to 7 angles using static light scattering at 35°, 50°, 75°, 90°, 105°, 130° and 145°. The detection takes place in a small cylindrical and high sensitivity light scattering cell with a total volume of 100 μ L. Due to the applied multi-angle detection, parameters such as absolute molar mass, radius of gyration and conformation can be determined. In this way, complex polymer and nanoparticle systems, where the use of single/dual angle light scattering technology is not sufficient, can be completely characterized. The PN3070 employs the latest opto-electronic components and devices, as e.g. modern microcontrollers, fiber optic detection coupled to a CCD and 24 bit A/D converters.

Unique Working Principle and Features

Avoiding the pitfalls of other commonly used MALS systems, the PN3070 brings the latest technology in electronics, CCD detection and light scattering cell design, directly into the scientific laboratory. Finally, there is a viable alternative available which offers unique advantages:

- Small footprint to save bench space and for more flexible instrument handling.
- Cylindrical cell design with vertical orientation to avoid sticking of airbubbles.
- Unique patented fiber optic detection without the use of fragile cell windows.
- Overcoming the index matching problem by using fiberoptic detection inside cell.
- Minimized dead volumes and absence of sharp narrow edges inside flow path.
- Extremely high dynamic range CCD detection overcoming outdated photodiodes.
- Direct gain adjustment via PC, no need for manual change of dip switches.
- The PN3070 is the new industry benchmark offered at a highly competitive pricing.



The patented fiberoptic flow cell incorporates 7 fused silica fibers, which are mounted directly in the cell showing only a low numerical aperture. There is no index matching necessary with this set-up and changes of the solvent does not effect the angular position of the measured stray light outside the cell, as it typically happens with other flow cell designs using conventional windows. The whole cell is designed for maximum ruggedness and for long-lasting reliable use, with pressures applied up to 35 bar.

Because of its simple and robust design, the cell and the instrument in total is ideally suited not only for R&D but also for QC applications monitoring production processes.

Typical Application Range

- Biotechnology/Pharma: Biopolymers, Starches, Proteins, Virus, Antibodies, Micelles
- Polymer Science: Plastics, Rubbers, Co-polymers, Latex particles, Polyelectrolytes
- Environmental Research: Organic and Inorganic Colloids, Humic and Fulvic Acids