**Ventillose
Taumelkolbenpumpe****REGLO-CPF Analog
Antrieb ISM1014****Kopf 2–25 µl/Hub**

FMI009	Keramikkolben
FMI010	Stahlkolben

Kopf 5–50 µl/Hub

FMI005	Keramikkolben
FMI013	LF = Low Flow

Kopf 10–100 µl/Hub

FMI007	Keramikkolben
--------	---------------

**Valveless Rotary
Piston Pump****REGLO-CPF Analog
Drive ISM1014****Head 2–25 µl/stroke**

FMI009	ceramic piston
FMI010	steel piston

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Piston pump-heads require thorough flushing after use

We recommend you flush the pump-head **after every pumping process** to prevent the piston from seizing by remaining particles (especially after pumping media containing salt, protein or particles).

Do not pump fluids containing suspended solids that are abrasive.

Please observe our note regarding filtration in the chapters Safety Precautions and Maintenance.

Piston fractures are not covered by our warranty.

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Caution

Rotating pumphead creates a pinch and crush hazard.

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Contents

Safety precautions	5
Warranty terms	9
Product	10
Rear panel	12
Mains voltage	13
Voltage setting	13
Changing the fuses	14
Starting the pump	14
Zero point calibration	15
Setting the stroke volume	16
Overload protector	17
Pumping against pressure	18
When the pump is not in use	18

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Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® gear pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- ▶ The pump must not be operated outside the designed operating and environmental conditions.
- ▶ The pump must not be used:
 - for medical applications
 - in explosion proof chambers or in the presence of flammable gases or fumes
- ▶ The circuit between mains supply and pump has to be connected to earth ground.
- ▶ Do not manipulate the pump-head before the pump is switched OFF and disconnected from the mains supply.
- ▶ The pump must be switched OFF before taking off the pump-head.

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Safety precautions

- ▶ Only new fuses, according to the specifications stated on Page 12 in this manual, must be used.
- ▶ The fuse-holder must not be short-circuited.
- ▶ Do not open or remove the housing while the pump is operating.
- ▶ Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken. When pumping corrosive media we recommend to place the pump in a collecting basin.
- ▶ Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- ▶ For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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Safety precautions

Media containing particles:

It is recommended that the pump-head be thoroughly flushed **after every pumping process** to prevent the piston from seizing due to remaining particles (especially after pumping media containing salt, protein or particles).

Do not pump media containing particles that exceed a diameter of 0.8 mm.

Piston stroke or stroke volume must be set large enough with respect to the particle size!

Piston fractures are not covered by our warranty.

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Safety precautions

Corrosive liquids:

Corrosive liquids may eventually produce leak paths around the sealing surfaces of the pump-head.

The pump should be inspected for leaks on a regular basis.

Dry running:

Before starting to run the pump it is recommended that the pump-head be filled with liquid.

Ensure that during priming the pump does not run dry for more than 15 sec.

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Warranty

For all parts manufactured by ISMATEC®

from date of delivery: 2 years
All other parts, excluding consumables: 1 year

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A pump must be returned in the original ISMATEC® packaging or in packaging of equal quality. The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded. Shipping costs are charged to the customer.

Our warranty becomes invalid in the case of:

- ▶ improper operation by the user, or if the pump is diverted from its proper use
- ▶ unauthorized modification or misuse by the user or by a third-party
- ▶ improper site preparation and maintenance
- ▶ operation outside of the environmental and electrical specifications for the product
- ▶ use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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Product

The REGLO-CPF Analog is a valveless rotary piston pump.

Package:

- ▶ Pump drive ISM1014
- ▶ Power Cord
- ▶ Operating Manual
- ▶ Pump-head as ordered

Head 2–25 µl/stroke

FMI009	ceramic piston
FMI010	steel piston

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Product

Please check the packaging and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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Rear panel

1. Analog interface
2. Mains socket
3. Fuse-holder with voltage selector 115/230V

Mains voltage

Mains voltage	Fuse rating
220–240 V _{ac}	2 x T800mA/250V
100–120 V _{ac}	2 x T800mA/250V

Socket/Power cord

Use exclusively the originally supplied power cord. The socket must be connected to the earth ground (protective conductor contact).

Before starting-up

Check that the voltage setting visible in the window of the fuse-holder complies with the local mains voltage. If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

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
Voltage setting 115V / 230V and changing the fuses

- ➔ Switch the pump OFF, pull out the
mains plug.
1. Pull out the fuse-holder by opening it with
a small screw-driver (size 0).
 2. Take out the voltage selector plate. Turn it
and re-insert it into the fuse-holder so that
the required voltage rating is acing the
window of the fuse-holder.
 3. Insert 2 new fuses
230 V_{ac}: 2 x T800mA/250V
115 V_{ac}: 2 x T800mA/250V
⚠ Always use 2 slow-blow fuses of the
same type complying with the local
mains voltage.
 4. Shut the fuse-holder. The voltage
rating is visible in the window.

⚠ Before starting-up

Check that the voltage setting visible in
the window of the fuse-holder complies
with the local mains voltage.

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
Starting the pump

1. Set the required speed on the 2-digit speed selector.
99 = 1800 rpm
(max. revolution)
➔ The rotation speed can be adjusted while the pump is running.
2. Choose the rotation direction.
3. Start the pump with the main switch = Run/Stop function.
4. Overload indicator (red LED).

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Zero point calibration

Before using the pump the first time, or whenever necessary, it is recommended that the 0-point of the flow control ring be calibrated.


- ▶ Start the pump and pump the liquid to be delivered.
- ▶ While the pump is in operation, turn the knurled adjustment ring to the right until the pump stops delivering the liquid (the motor is still running!).
- ▶ Loosen the fixing screw of the revolving flow control ring (black) and turn it to the value '0'.
- ▶ Tighten the fixing screw again.

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Setting the stroke volume

For setting the stroke volume, turn the knurled adjustment ring till the desired value (0–450) on the fixed and revolving flow control ring appears. (The fixing screw must not be loosened!)

Basically the following applies:

- ▶ **Small stroke volume/ high speed:**
 - low pulsation
 - better resolution of the flow rate and dispensing volume
 - tendency to cavitation (insufficient flow-through of the liquid in the pump-head).
- ▶ **Large stroke volume/ low speed:**
 - for elevated viscosity
 - for differential pressures
 - reduced resolution of the flow rate.

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Overcurrent protector

The pump REGLO-CPF Analog features an overload protector. When an overload condition occurs, an activated overcurrent protector is indicated by the red LED and the pump is stopped.

Whenever this situation occurs, **the pump must be switched OFF immediately.**

➔ Let the pump cool down (it takes about 2 minutes until the pump is ready again).

Before the pump is re-started, it is important to check the reason for the overload (eg. too high differential pressure).

➔ Only after the cause of the overload has been detected and the failure corrected accordingly, may the pump be started again.

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Pumping against pressure

The REGLO-CPF Analog can be used for continuous duty at a differential pressure of max. 6.9 bar (100 psi).

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
When the pump is idle

It is recommend that the pump-head be flushed **after every pumping process** to prevent the piston from seizing by remaining particles (especially after pumping media containing salt, protein or particles).

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Pin 2, remote
Pin 3, start
Pin 4, direction
Pin 13, internal speed

Analog-Eingang
Analog input/ Entrée analogique
Pin 5, speed IN
0-5 VDC / 0-10 VDC / 0-20 mA / 4-20 mA

Analog-Ausgang (TTL-Pegel)
Analog output (TTL level)
Sortie analogique (niveau TTL)
Pin 9, speed OUT 0-9 kHz

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Analog interface

Pin 1, GND (ground)
Reference potential for all other inputs.

Pin 2, remote
For changing between manual control and analog interface. For activating the analog interface, Pin 2 must be connected with Pin 1 (GND).

Pin 3, start
In remote operation (Pin 2 to GND) the pump starts when connected to Pin 1 (GND).

Pin 4, direction
In the open position the pump turns counter-clockwise; when connected to Pin 1 (GND) it turns clockwise.

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
Digitale Eingänge (TTL-Pegel)

Digital inputs (TTL-level)

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- Pin 2, remote
- Pin 3, start
- Pin 4, direction
- Pin 13, internal speed

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Analog interface

Pin 5, speed IN

For external speed control (0–5V_{DCr}, 0–10V_{DCr}, 0–20mA, 4–20mA) Input impedance and alternatives by means of DIP switches inside the pump (see page 22).

Pin 7, +22V_{DC}

About +22 V_{DC} are available (max. current 0.2 A).

Pin 9, speed OUT

0–10 kHz Frequency proportional to the rotation speed: 0–5000 rpm (with each revolution of the piston, 300 impulses are generated).

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
Analog interface

Pin 10, +5V_{cc}

About +5 V_{cc} are available.
 (max. current 0.1 A).

Pin 13, internal speed

Analog interface activated
 (Pin 2 on GND).

- ▶ Pin 13 open: The rotation speed is adjusted via Pin 5 (speed IN).
- ▶ Pin 13 on GND: The rotation speed can be adjusted by the speed selector on the control panel of the pump.

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
Settings of switch S1

➔ Make sure that the pump is disconnected from the mains supply.

⚠ The instrument should only be opened by a qualified technician!

Capacitors inside the pump may still be charged even though the mains plug has been disconnected for some time.

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Pins		Impedanz	DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 5 speed IN	0-5V _{dc}	18kΩ	ON*	OFF*	OFF*	OFF*	OFF*	ON*
	0-10V _{dc}	38kΩ	OFF	ON	OFF	OFF	OFF	ON
	0-20mA	250Ω	OFF	OFF	ON	OFF	ON	ON
	4-20mA	250Ω	OFF	OFF	OFF	ON	ON	OFF

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Footswitch

Order No. ISM891

This footswitch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc. Both hands are free for handling the bottles and tubing.

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Please note

In order to reach the max. flow rate, the stroke volume must be set to the max. value (25, 50 or 100 µl) depending on the pump-head model.

If, depending on the application, this stroke volume must be reduced, it is not recommended the adjustment go below 10% of the maximum value.

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*Range is outside the recommended min. stroke volume

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Maintenance

It is recommended that the the pump-head be thoroughly flushed after pumping to prevent the piston from seizing due to remaining particles (especially when media has been pumped that contained salt, proteins or particles).

Otherwise, the REGLO-CPF Analog requires no particular maintenance or mechanical adjustments.

Avoid pumping media containing abrasive solids.

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Disposal

Please retain packaging materials until the product warranty ends. Afterwards please discard packaging materials in an environmentally-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Service and repairs

Like all pumps, the REGLO-CPF Analog contains components that will wear over a period of time.

For repairs, send the complete pump with detailed description of the failure to your ISMATEC® distributor. Please use the original ISMATEC® packaging or a packaging of equal quality.

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Please Note

For maintenance, service, repairs and warranty of the pump-head we also refer to the enclosed instruction manual of the pump-head manufacturer. In case of doubt that manual shall prevail.

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Replacement parts

For repairs performed by the owner (out of the warranty period), ask your ISMATEC® distributor for:

► Piston-/cylinder group

Type	Order No
H-00CKC-LF	FMI311
H-00SKY-LF	FMI310
H-0CKC	FMI307
H-1CKC	FMI308

► Lip seals (2 pieces are needed)

Type	Order No
H408-00J (RH00)	FMI325
H408-0A (RH0)	FMI312
H408-1A (RH1)	FMI319

- Parts lists
- Replacement parts for pump drive
- Wiring diagrams

Please give information on defect, date of purchase, serial-no., model.

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Speed

18–1800 rpm, digitally adjustable in steps of 1%.

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Please Note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Technical Specifications

Operating conditions

Temperature +5°C to +40°C
(+41°F to 104°F)
Rel. humidity max. 80%
not condensing, at normal
laboratory conditions

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
Operating conditions

1. REGLO-CPF drive
2. Coupling completely assembled
3. Mounting plate
4. Countersunk head-screws (4 pcs.)
5. Allen Key 2.5 mm

from FMI

- a. CPF pump-head (RH00, RH0, RH1)
- b. Screws (2 pieces)
- c. Serrated lock washers (2 pieces)

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


Mounting a ceramic piston pump-head on the REGLO-CPF drive

1. Attaching the mounting plate

Remove screws (b) and serrated lock washers (c) from the pump-head.

- Fasten the mounting-plate (3) on the pump-head with the 2 screws (b) and serrated lock washers (c).

Important

The 4 countersinks (5) must be concave to the pump-head.

3. Fixing the coupling

Place the completely assembled coupling (2) onto the pump-head shaft and push it on as far as it will go.

The Allen screw (d) must be pointing towards the flat part of the pump-head shaft (e).

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4. Tighten the coupling (2) with the Allen Key size 2.5 mm (5).

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
Important

- The scale of the red flow control ring (f) must be pointing upwards.
- Hold the mounting plate in such a way that the screws (b) are on the top and at the bottom.
- Turn the coupling so that afterwards the Allen screw (d) can be tightened through the aperture (g) with the Allen Key (5).
- Insert the pump-head assembly and push it on as far as it will go. Check the position of the Allen screw at the aperture (g).


1. Drive in the flat head machine-screws (4). Tighten the Allen screw (d) with an Allen Key (5) size 2.5 mm.

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6

ISMATEC®

Learn more
about the
Ismatec® product
line by visiting:
www.ismatec.com

For ordering and technical support, please contact:

North America sales@ismatec.com | 1-800-323-4340 | 1-847-549-7600

Europe sales.europe@ismatec.com | +49 (0) 9377 9203-0

 **Cole-Parmer®**