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Gear pump drive with
variable speed drive

BVP-Z Standard

Drive without pump-head
Order No. ISM446

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Caution

Rotating pumphead creates a pinch and crush hazard.

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Contents

Safety precautions	3 – 4
Warranty terms	5
Product	6
Operating panel	7
Rear panel	7
Mains voltage	8
Voltage setting	8
Changing the fuses	8
Starting the pump	9
Overcurrent protector	10
Analog interface	11 – 13
Pumping against pressure	14
Pump-head with bypass	14
Before switching the pump off	14
Pump-head	15
Mounting the pump-head	15 – 16
Gear Pump Technology	17
Pump-head specifications	17
Pump-heads, overview	18 – 19
Accessories	
Centering flange	20
Non-return valve	20
Tubing clamps	20
Tubing nipple	20
Footswitch	21
Service Kit	21
Maintenance, repairs	22
Spare parts	22
Disposal	22
Technical specifications	23
Mounting instructions	24
Centering flange for pump-heads Micropump® 200 and 201	

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Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® gear pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and the following recommendations will be observed.

- ▶ The circuit between mains supply and pump has to be earthed.
- ▶ The pump must not be operated outside the destined operating and environmental conditions.
- ▶ Place the unit in a well ventilated position at least 10 cm away from walls, partitions etc. Ensure that curtains and similar materials do not cover the ventilation slits.
- ▶ The pump **must not** be used:
 - for medical applications on human beings
 - in explosion-proof chambers or in the presence of flammable gases or fumes.
- ▶ The pump must be switched OFF when changing the pump-head.
- ▶ The permeability of pump tubing depends on the material used and pressure conditions. Tubing can also become electrostatically charged. Please be aware of possible hazards when routing tubing in explosion-proof chambers.

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Please note

Never pump media containing particles.

If necessary, filter the media before starting your application.

Never run the pump-head under dry conditions.

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Safety precautions

- ▶ Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- ▶ Do not manipulate the pump-head before the pump is switched off and disconnected from the mains.
- ▶ Before initial start-up, the pump-head must be filled with the pumping liquid or water. Never run the pump-head under dry conditions.
- ▶ We advise you never pump liquid containing particles. Even the smallest particles have an abrasive effect, which can ruin the gears in the course of a few seconds.
- ▶ Only new fuses, according to the specifications stated on Page 8 in this manual, must be used.
- ▶ The fuse holder must not be short-circuited.
- ▶ Do not open or remove the housing while the pump is operating.
- ▶ Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- ▶ For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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BVP-Z drive: 1 year

Pump-head: 1 year

In case of any queries, please contact your nearest ISMATEC® representative.

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Warranty Terms

We warrant the perfect functioning of the BVP-Z drive, provided it has been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

The pump-heads are covered by the specific terms of Micropump®. Gears and seals are regarded as consumables and as such they are not covered by the warranty.

The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded. Shipping costs for instruments under warranty are charged to the customer.

Our warranty becomes invalid in the case of:

- ▶ improper operation by the user, or if the pump is diverted from its proper use,
- ▶ unauthorized modification or misuse by the user or by a third-party,
- ▶ improper site preparation and maintenance,
- ▶ operation outside environmental and electrical specifications for the product,
- ▶ use of third-party software, hardware, interfacing or consumables purchased by the user and which do not comply with our specifications,
- ▶ damages caused by contamination or leaks due to torn or burst tubing.

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Product

Contents of the package

- ▶ BVP-Z Standard drive
1 magnet
Order No. ISM446
- ▶ 2 tubing clamps
- ▶ 1 non-return valve (PP)
- ▶ Operating manual

If ordered:

- ▶ Pump-head with
2 tubing adapters
- ▶ other accessories

Please check the package

and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

➡ Complaints can only be accepted within 8 days from receipt of the goods.

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

Operating panel

1. Mains switch
2. Control light
3. Speed selector
999 = 6000 rpm (max. speed)
= 100 rps adjustable in 0.1% steps
4. RUN/STOP
5. MAX key for drive speed (ideal for fast filling or emptying the system)

⚠ **The drive runs only in clock-wise direction.**

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
Mains voltage

Mains voltage	Voltage setting	Fuse rating
220–240 V _{AC} 110–120 V _{CA}	230 V 50/60Hz 115 V 50/60Hz	2 x T1.25A/250V 2 x T2.50A/250V

⚠ Before starting up

Check if the voltage setting visible in the window of the fuse holder complies with your local mains voltage.

If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

⚠ Socket/Power cord

Use exclusively the original power cord supplied with the pump. The socket must be earthed (protective conductor contact).

Voltage setting 115V_{AC} / 230V_{AC} and changing the fuses

➡ Switch the pump off,
pull out the mains plug.

1. Pull out the fuse holder by opening it with a small screw-driver (size 0).
2. Take out the voltage selector plate. Turn it and re-insert it into the fuse holder so that the required voltage rating is facing the window of the fuse holder.
3. Insert 2 new fuses
➡ 230 V: 2 x 1.25 A (slow-blow)
➡ 115 V: 2 x 2.50 A (slow-blow)

⚠ Use always 2 slow-blow fuses of the same type complying with the local mains voltage. In case of a defective fuse, always replace both.

4. Shut the fuse holder. The voltage rating is visible in the window.

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Please note


Before initial start-up, the pump-head must be filled with the medium or water. Never run the pump-head under dry conditions. Never pump liquids that contain particles. If necessary, filter the medium before use.

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Starting the pump

- ▶ Check the voltage rating in the window of the fuse holder (on rear panel). If necessary, change the rating as indicated on Page 8.
- ▶ If necessary, mount the pump-head as stated on Page 15.
- ▶ Before screwing the tubing nipples into the pump-head, we recommend you wrap PTFE tape around the thread of the nipple. This ensures optimum sealing and helps to protect the thread.
- ▶ Mount the tubing and connect the hose line to your system. Fasten the tubing to the nipples by means of the tubing clamps.
- ▶ Connect the pump to the mains.
- ▶ Fill the pump-head with liquid before you start pumping.
 1. Press in the mains switch
 2. The control light is on
 3. Set the required speed on the 3-digit speed selector. (999 = 6000 rpm)
 - ➔ The rotation speed can be adjusted while the pump is running.
 4. Start the pump by switching the RUN/STOP button.
 5. The MAX button allows rapid filling (or emptying) of the system while the pump is running.

⚠ **Before you start pumping, fill the pump-head with the medium or water. Never run the pump-head under dry conditions.**

⚠ **Never pump liquids that contain particles. If necessary, filter the medium before use.**

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The tubing pump BVP *Standard* features an overload protector.

When the overcurrent protector is activated, the pump is stopped and the control light flashes (see the description of Pin 12 of analog interface on Page 12).

Whenever this situation occurs, **the pump must be switched OFF immediately.**

➡ Let the pump cool down (it takes about 2 minutes until the pump is ready again).

Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure).

➡ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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Analog Interface

Pin 9, speed OUT

The default setting is 0–10 V_{DC}, proportionally to the motor speed 0–6000 rpm. Alternatively a frequency range from 0/12 kHz is available. The input range can be selected via a DIP-switch inside the pump (see Page 13).

Pin 10, +5V_{DC}

This pin is connected to the stabilized +5V_{DC} voltage supply of the pump.

Pin 12, idle, overload

This pin is at +5V_{DC} while the pump runs. When idle or in overload, the voltage changes to 0V_{DC}.

➔ Switch off the mains voltage and let the device cool down for 2 minutes.

Pin 13, internal speed

If Pin 13 is connected with Pin 1 (ground), the pump speed can be adjusted via the speed selector on the operating panel instead of using the signal at Pin 5.

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
- A = 0 – 10 V_{DC} (Standard)
- B = 0/12 kHz

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Settings of switch S1

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Pins			DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 5 speed IN	0-5V	100kΩ	ON*	OFF*	OFF*	OFF*	OFF*	ON*
	0-10V	20kΩ	OFF	ON	OFF	OFF	OFF	ON
	0-20mA	250Ω	OFF	OFF	ON	OFF	ON	ON
	4-20mA	250Ω	OFF	OFF	OFF	ON	ON	OFF

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➡ **Make sure that the pump is disconnected from the mains supply.**

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⚠ **The instrument should only be opened by a qualified technician!**
Capacitors inside the pump may still be charged even though the mains plug has been disconnected for some time.

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Pump-head

In general, the BVP-Z Standard gear pump is supplied as a **complete unit**, consisting of: the BVP-Z Standard drive and

1. driving magnet and as ordered:
2. pump-head Micropump®
3. tubing adapters (2 pcs.)

For the BVP-Z Standard, a wide range of different gear pump-heads with PTFE, PPS or graphite gears are available.

You find an overview chart of the pump-heads available on Pages 18/19.

⚠ **Before mounting a new pump-head, make sure that the correct magnet is fitted to the drive.**

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Mounting the pump-head

Fig. 2

Fasten the driving magnet on the motor shaft (through the hole underneath).

➔ Ensure that there is a gap of about 3 mm between the magnet and the brim of the fixing case.

⚠ **During operation the back of the magnet must not touch the housing of the pump drive.**

Fig. 3

Screw the 2 tubing adapters into the pump-head.

➔ Before screwing the tubing adapters into the pump-head, it is recommended that PTFE tape be wrapped around the thread of the adapter. This ensures optimum sealing and helps to protect the thread.

Fig. 4

Screw the pump-head to the fixing case of the drive. The face of the pump-head with the engraved model number must be exposed to the user. The engraved arrow indicates the flow direction.

➔ Please note that gear pumps can only pump in one direction. An arrow engraved on the front of the pump-head indicates the correct flow direction.

⚠ **Fill the pump-head with liquid before you start pumping. Avoid dry-running of the pump-head.**

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Fig. 2

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»Cavity style« pump-head 

- ▶ The »Cavity style« pump-head can be used for viscous media and for applications with a certain suction lift
- ▶ Excellent chemical resistance
- ▶ Smooth operation at a low noise
- ▶ Low inner friction
- ▶ Incl. B mount (flange)

»Suction shoe« pump-heads 

- ▶ Enhanced pump performance at elevated differential pressures
- ▶ Fewer parts, easy to service
- ▶ Suitable for elevated temperature range
- ▶ Not suitable for applications requiring a suction lift
- ▶ Incl. B mount (flange)

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Pump-head specifications

- ▶ Stainless steel pump-head housing SS 316/DIN 1.4571

▶ **Temperature range**

- 46°C to +54°C / -51°F to +129°F (Cavity style)
- 46°C to +177°C / -51°F to +350°F (Suction shoe)(higher temperatures possible on request)

▶ **Differential pressure**

- The difference between the pressure on the inlet and the discharge port of the pump-head max. 5.2 bar (75 psi) depending on pump-head (Pages 18/19)

▶ **System pressure**

- The pressure within the pump-head max. 21 bar / 305 psi

▶ **Gear material**

- PTFE = Polytetrafluoroethylene
- PPS = Polyphenylenesulphide
- Graphite = Carbon graphite
- PEEK = Polyetheretherketone

▶ **Ports**

- (Inside thread) 1/8"-27 NPT

The pump-heads are interchangeable. Different drives take different magnets. The BVP-Z Standard drive is originally supplied with the correct magnet.

When ordering a new pump-head, please state the intended operation temperature. If you intend to use a Cavity style pump-head for temperatures exceeding 54°C (130°F), the correct sealing material must be ordered. Please contact your nearest ISMATEC® representative. Some pump-heads are available with an **internal bypass**. This bypass is to be used for safety reasons only and must not be used as pressure control (see Page 14).

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Footswitch

Order No. ISM891

This footswitch serves as a start/ stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc. Both hands are free for handling the bottles and tubing.

Service Kit

The service kit provides all spare parts in case that gears and seals require replacement.

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Disposal

Please retain packaging materials until the product warranty ends. Afterwards please discard packaging materials in an environmentally friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Maintenance

After use, it is recommended that the pump-head be thoroughly flushed with warm tap water.

Provided the BVP-Z *Standard* is operated properly and in compliance with this manual, the moving parts such as gears, etc. are the only parts that are subject to wear and tear.

The Service Kit corresponding to the individual pump-head allows the user to replace the parts that are subject to wear and tear (see Page 21).

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Repairs

For repairs please send the BVP-Z *Standard* drive to your ISMATEC® representative.

Please give information on defect, date of purchase, Serial No., and model.

To request a repair, please fill out the form you can download here: www.ismatec.com

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Spare parts

For repairs carried out by the owner (out of the warranty period) please ask your ISMATEC® distributor for:

- ▶ spare parts
- ▶ parts lists
- ▶ wiring diagrams

Please give information on defect, date of purchase, Serial No., and model.

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Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Max. 150 W

Operating conditions

Temperature +5 to +40°C,
(+41 to 104°F) Rel.
humidity max. 80%
not condensing, at normal
laboratory conditions

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

Centering Flange

The Centering Flange fits onto the BVP-Z drive. It is required only for the Micropump® 200 and 201 pump-head series where it facilitates their precise mounting.

Accessories:

- 4 flat-headed screws
M 2.5 x 8 BN0187
- 1 centering flange
 IS10379

Mounting instructions:

- Please ensure that the circular drive magnet (c) is correctly mounted.
 - When mounting the centering flange, make sure that
 - its elevated inner ring (1) exactly fits into the drive flange, and
 - the two straight holes (2) are placed in a vertical position at the top and the bottom
- ➡ When screwing the centering flange onto the drive, make certain that none of the four flat-headed screws protrude.
- You can now proceed with mounting the gear pump-head with the two screws in the usual way.

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ISMATEC®

Learn more
about the
Ismatec® product
line by visiting:
www.ismatec.com

For ordering and technical support, please contact:

North America sales@ismatec.com | 1-800-323-4340 | 1-847-549-7600

Europe sales.europe@ismatec.com | +49 (0) 9377 9203-0

 **Cole-Parmer®**