

# EAF2000 Electrical Flow FFF Series

Electrical Asymmetrical Flow FFF




**Universal Particle/Polymer/Protein Separator**

# EAF2000 Electrical Flow FFF Series

## Simultaneous Electrical and Asymmetrical Flow FFF

### Principle of EAF4

The EAF4 technology combines the principle of Electrical and Asymmetrical Flow FFF in one system. Electrical and Cross Flow fields are applied simultaneously across the FFF channel. Separations by particle size and particle charge based on electrophoretic mobility can be achieved.



EAF2000 Electrical Flow FFF Channel



PN2410 Electrical FFF Module



### Hardware & Software

The EAF2000 System is based on the AF2000 Asymmetrical Flow FFF module, with additional Electrical FFF module, special EAF4 channel and software module:

- PN2410 Electrical FFF module for precise control of the Electric Field allowing constant and reproducible separation conditions
- EAF2000 Electrical Flow FFF channel with special built-in pole-reversible electrodes and resistant membrane for cross flow
- Software module for constant current operation and data evaluation for electrophoretic mobility calculation

## Specifications

### EAF4 Unit

- Carrier Liquids:
  - Aqueous: any aqueous liquid, pH from 2 - 11, ionic strength from DI water to saline
  - Organic: THF, MeOH, etc.
- Detectors: UV, RI, laser light scattering, viscometry, fluorescence and more; multiple detectors in-line
- System Software: NovaFFF AF2000 Control, NovaFFF Analysis
- Temperature Range: Room temperature
- PC Requirements: Windows, min. 8 GB RAM, 2 Ethernet-LAN ports, 1 RS232 port

### PN2410 Electrical FFF Module

- Output / E-Field:
  - Current: +/- 75 mA \*
  - Max. voltage: +/- 22.5 V
  - Resolution: 10  $\mu$ A
  - Precision: <1 % \*\*
  - Reproducibility: <0.5 % RSD \*\*
- Status Control:
  - Current read back
  - Voltage read back
  - Polarity
  - Discharge
- User Control:
  - Current
  - Selection of polarity
  - Electrode discharge to make sure E-field is zero
- Safety:
  - Short-circuit proof
- Status Indications LED:
  - Red = field OFF
  - Green = field active
- Environmental Conditions:
  - Relative humidity 20 – 80 % (non-condensing) at an operating temperature range of 10 – 40 °C
- Power Requirements:
  - 100 – 240 VAC @ 50 – 60 Hz
- Dimensions (DxWxH):
  - 430 x 270 x 90 mm
- Weight:
  - 5 kg

\* Depends on used solvents and channels. Very low conductivities will cause lower currents.  
 \*\* Valid only for an output current between 0.5 mA and 70 mA.

### Applications of EAF2000

Biopharm : Proteins, Antibodies, Antibody Aggregates, Viruses, Liposomes and Drug Delivery  
 Environment : Environmental Nanoparticles such as Humics, Fulvics, Clay Particles, Silica, TiO<sub>2</sub>  
 Nanomaterials : Polystyrene Latex Beads, any charged Nanoparticle or high-tech Nanomaterial

The Postnova EAF2000 Electrical Flow FFF technology is compatible with any existing AF2000.

Technical specifications are subject to change without further notice.

## Contact

- Postnova Analytics GmbH  
86899 Landsberg, GERMANY  
T: +49 8191 985 688 0
- Postnova Analytics UK Ltd.  
Malvern, Worcestershire, WR14 3SZ, UK  
T: +44 1684 585167
- Postnova Analytics Inc.  
Salt Lake City, UT 84102, USA  
T: +1 801 521 2004
- Postnova North Europe  
01630 Vantaa, FINLAND  
T: +358 9 8545 510

info@postnova.com  
 www.postnova.com

