

POLYMERSTANDARDS




Table of Contents

Polymer Kits

Organic Kits	4
Aqueous Kits	4
DIN Kits	4
ReadyCal Kits	5
Validation Kits	5
MALDI Kits	5

Organic Polymers

Nylon 6 broad standard	6
Nylon 66 broad standard	6
Oligo(ethylene) monodisperse	6
Poly(2-vinyl pyridine)	6
Poly(acryl nitrile) broad standard	7
Poly(alpha-methyl styrene)	7
Poly(butadiene 1,2)	7
Poly(butadiene 1,4)	8
Poly(butadiene d6)	8
Poly(butylene terephthalate) broad standard	8
Poly(carbonate)	9
Poly(carbonate)-d4 deuterated	9
Poly(dimethyl siloxane)	9
Poly(dimethyl siloxane)-d6 deuterated	10
Poly(ethyl methacrylate)	10
Poly(ethylene terephthalate) broad standard	10
Poly(ethylene) broad standard	10
Poly(iso-butylene)	11
Poly(iso-propyl methacrylate)	11
Poly(isoprene 1,4)	12
Poly(isoprene 3,4)	12
Poly(methyl methacrylate)	13
Poly(methyl methacrylate) broad standard	14
Poly(methyl methacrylate) isotactic	14
Poly(methyl methacrylate) syndiotactic	14
Poly(methyl methacrylate)-d8 deuterated	14
Poly(n-butyl methacrylate)	15
Poly(n-propyl methacrylate)	15
Poly(o-chloro styrene)	15
Poly(p-methyl styrene)	15
Poly(p-methyl styrene) deuterated	15
Poly(propylene) broad standard	15
Poly(styrene)	16
Poly(styrene) DIN cert.	17
Poly(styrene) broad standard	17
Poly(styrene)-d8 deuterated	18

Poly(t-butyl acrylate)	18
Poly(tert.butyl methacrylate)	19
Poly(tetrahydrofuran)	19
Poly(vinyl acetate) broad standard	20
Poly(vinyl butyral) broad standard	20
Poly(vinyl chloride)	20
Polyethylene	20
Polypropylene	21
Aqueous Polymers	
Dextrane	22
Dextrane DIN cert.	22
Dextrane broad standard	22
Hydroxy-ethyl starch	22
Hydroxy-propyl cellulose	23
Hydroxy-propylcellulose broad standard	23
Poly(2-vinyl pyridine) quaternated	23
Poly(acrylamid) broad standard	24
Poly(acrylic acid) sodium salt	24
Poly(DADMAC)	24
Poly(allyl amine) hydrochloride	24
Poly(ethylene glycol)	25
Poly(ethylene glycol) DIN cert.	25
Poly(ethylene oxide)	26
Poly(ethylene oxide)-d4 deuterated	26
Poly(methacrylic acid) sodium salt	26
Poly(methacrylic acid)	26
Poly(propylene glycol)	27
Poly(styrene sulfonate) sodium salt	27
Poly(vinyl alcohol) broad standard	27
Poly(vinyl pyrrolidone) broad standard	28
Pullulane	28

Polymer Kits

Organic Kits

Article No.	Description	Mp [D]	Amount
Z-POS-nsskit	Nylon 66 broad standard	22.000-70.000	7 x 0,25 g
Z-POS-oekit	Oligo(ethylene) monodisperse	112-2.200	8 x 0,50 g
Z-POS-pvpkit	Poly(2-vinyl pyridine)	1.000-300.000	8 x 1 g
Z-POS-amskit	Poly(alpha-methyl styrene)	1.500-1.000.000	10 x 1 g
Z-POS-bdkit	Poly(butadiene 1,4)	1.000-1.000.000	10 x 1 g
Z-POS-pckit	Poly(carbonate)	500-180.000	8 x 0,5 g
Z-POS-pdmkit	Poly(dimethyl siloxane)	400-300.000	8 x 0,5 g
Z-POS-emkit	Poly(ethyl methacrylate)	2.000-500.000	10 x 1 g
Z-POS-pibkit	Poly(iso-butylene)	100-1.000.000	10 x 0,5 g
Z-POS-piokit	Poly(isoprene 1,4)	1.000-800.000	10 x 1 g
Z-POS-mmkitl	Poly(methyl methacrylate)	100-30.000	6 x 1 g
Z-POS-mmkith	Poly(methyl methacrylate)	100-1.000.000	12 x 1 g
Z-POS-nbkit	Poly(n-butyl methacrylate)	1.500-800.000	9 x 1 g
Z-POS-pskitl	Poly(styrene)	100-20.000	8 x 1 g
Z-POS-pskith	Poly(styrene)	400-3.000.000	12 x 1 g
Z-POS-tbkit	Poly(tert.butylmethacrylate)	600-1.000.000	10 x 1 g
Z-POS-ptkith	Poly(tetrahydrofurane)	200-100.000	8 x 1 g
Z-POS-pvckit	Poly(vinyl chloride)	100-320.000	10 x 0,5 g
Z-POS-pekit	Polyethylene	100-170.000	10 x 0,25 g

Aqueous Kits

Article No.	Description	Mp [D]	Amount
Z-POS-dxtkit	Dextrane	180-300.000	10 x 0,5 g
Z-POS-heskit	Hydroxy-ethyl starch	9.000-2.000.000	10 x 0,25 g
Z-POS-paakit	Poly(acrylic acid) sodium salt	1.700-1.100.000	10 x 0,25 g
Z-POS-pegkit	Poly(ethylene glycol)	100-40.000	10 x 1 g
Z-POS-peokit	Poly(ethylene oxide)	20.000-1.700.000	8 x 0,5 g
Z-POS-pmakit	Poly(methacrylic acid) sodium salt	1.000-1.000.000	10 x 1 g
Z-POS-ppgkit	Poly(propylene glycol)	80-5.000	7 x 1 g
Z-POS-psskit	Poly(styrene sulfonate) sodium salt	4.000-1.000.000	10 x 1 g
Z-POS-pvokit	Poly(vinyl alcohol) broad standard	5.000-140.000	8 x 1 g
Z-POS-pulkitl	Pullulane	5.000-800.000	8 x 0,2 g
Z-POS-pulkith	Pullulane	300-800.000	10 x 0,2 g

DIN Kits

Article No.	Description	Mp [D]	Amount
Z-POS-dxtkitd	Dextrane	180-300.000	10 x 0,5 g
Z-POS-pegkitd	Poly(ethylene glycol)	100-40.000	10 x 1 g
Z-POS-mmkitd	Poly(methyl methacrylate)	100-1.000.000	12 x 1 g
Z-POS-pskitd	Poly(styrene)	200-2.000.000	12 x 1 g

ReadyCal Kits

Article No.	Description	Mp [D]	Amount
Z-POS-pegkitr1	PEG ReadyCal 1,5 mL Vials	100-40.000	3 x 10
Z-POS-peokit1	PEO ReadyCal 1,5 mL Vials	200-1.200.000	3 x 10
Z-POS-peokit4	PEO ReadyCal 4 mL Vials	200-1.200.000	3 x 10
Z-POS-mmkitr1	Poly(methyl methacrylate), ReadyCal 1,5 mL Vials	500-3.000.000	3 x 10
Z-POS-mmkitr4	Poly(methyl methacrylate), ReadyCal 4 mL Vials	500-3.000.000	3 x 10
Z-POS-pskitr1	Poly(styrene) ReadyCal 1,5 mL Vials	400-3.000.000	3 x 10
Z-POS-pskitr4	Poly(styrene) ReadyCal 4,0 mL Vials	400-3.000.000	3 x 10
Z-POS-pskitr1l	Poly(styrene) ReadyCal 1,5 mL Vials	162-70.000	3 x 10
Z-POS-pskitr4l	Poly(styrene) ReadyCal 4,0 mL Vials	162-70.000	3 x 10
Z-POS-pskitr1h	Poly(styrene) ReadyCal 1,5 mL Vials	1.500-6.500.000	3 x 10
Z-POS-pskitr4h	Poly(styrene) ReadyCal 4,0 mL Vials	1.500-6.500.000	3 x 10

Validation Kits

Article No.	Description	Mp [D]	Amount
Z-POS-dxkit	Dextran	10.000-230.000	4 x 0,5 g
Z-POS-mmkit	Poly(methylmethacrylat)	8.000-1.100.000	4 x 0,5 g
Z-POS-pskit	Poly(styrol)	9.000-800.000	4 x 0,5 g

MALDI Kits

Article No.	Description	Mp [D]	Amount
Z-POS-pegkitm	Poly(ethylene glycol)	500-20.000	6 x 0,5 g
Z-POS-mmkitm	Poly(methyl methacrylate)	500-70.000	6 x 0,5 g
Z-POS-pskitm	Poly(styrene)	500-70.000	6 x 0,5 g

Organic Polymers

Nylon 6 broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-ns17k	1	17.000	<1,80
Z-POS-ns20k	1	20.000	<1,80
Z-POS-ns30k	1	30.000	<1,80
Z-POS-ns40k	1	40.000	2,11

Nylon 66, broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-nss30k	0.25	30.000	<3,00
Z-POS-nss35k	0.25	35.000	<3,00
Z-POS-nss40k	0.25	40.000	<3,00
Z-POS-nss50k	0.25	50.000	<3,00
Z-POS-nss70k	0.25	70.000	<3,00
Z-POS-nss110k	0.25	110.000	<3,00

Oligo(ethylene) monodisperse

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-oe72	1	72	1,00
Z-POS-oe86	1	86	1,00
Z-POS-oe114	1	114	1,00
Z-POS-oe142	1	142	1,00
Z-POS-oe170	1	170	1,00
Z-POS-oe200	1	200	1,00
Z-POS-oe225	1	225	1,00
Z-POS-oe280	1	280	1,00
Z-POS-oe340	1	340	1,00
Z-POS-oe390	1	390	1,00
Z-POS-oe450	1	450	1,00
Z-POS-oe500	1	500	1,00
Z-POS-oe600	1	600	1,00

Poly(2-vinyl pyridine)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pvp1.2k	1	1.200	<1,10
Z-POS-pvp3k	1	3.000	1,18
Z-POS-pvp4.5k	1	4.500	<1,10
Z-POS-pvp6.5k	1	6.500	<1,10
Z-POS-pvp11k	1	11.000	<1,10
Z-POS-pvp12k	1	12.000	<1,10
Z-POS-pvp14k	1	14.000	<1,10
Z-POS-pvp20k	1	20.000	<1,10
Z-POS-pvp65k	1	65.000	<1,10
Z-POS-pvp80k	1	80.000	<1,10
Z-POS-pvp115k	1	115.000	<1,10
Z-POS-pvp135k	1	135.000	1,13
Z-POS-pvp260k	1	260.000	1,26
Z-POS-pvp290k	1	290.000	<1,10
Z-POS-pvp350k	1	350.000	<1,10
Z-POS-pvp850k	1	850.000	<1,10
Z-POS-pvp1.2m	1	1.200.000	1,37

Poly(acryl nitrile) broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pan30k	0.5	30.000	<3,00
Z-POS-pan35k	0.5	35.000	<3,00
Z-POS-pan65k	0.5	65.000	<3,00

Poly(alpha-methyl styrene)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-ams1.5k	1	1.500	1,29
Z-POS-ams2k	1	2.000	1,12
Z-POS-ams4k	1	4.000	<1,10
Z-POS-ams8k	1	8.000	<1,10
Z-POS-ams15k	1	15.000	1,12
Z-POS-ams35k	1	35.000	<1,10
Z-POS-ams40k	1	40.000	<1,10
Z-POS-ams60k	1	60.000	<1,10
Z-POS-ams70k	1	70.000	<1,10
Z-POS-ams100k	1	100.000	<1,10
Z-POS-ams110k	1	110.000	<1,10
Z-POS-ams200k	1	200.000	<1,10
Z-POS-ams250k	1	250.000	<1,10
Z-POS-ams400k	1	400.000	<1,10
Z-POS-ams900k	1	900.000	<1,10

Poly(butadiene 1,2)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-bdt800	1	800	<1,10
Z-POS-bdt1.2k	1	1.200	<1,10
Z-POS-bdt1.5k	1	1.500	<1,10
Z-POS-bdt2k	1	2.000	<1,10
Z-POS-bdt2.5k	1	2.500	<1,10
Z-POS-bdt3.5k	1	3.500	<1,10
Z-POS-bdt4k	1	4.000	<1,10
Z-POS-bdt6k	1	6.000	<1,10
Z-POS-bdt12.5k	1	12.500	<1,10
Z-POS-bdt100k	1	100.000	<1,10

Poly(butadiene 1,4)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-bdf110	1	181	<1,10
Z-POS-bdf350	1	350	<1,10
Z-POS-bdf700	1	700	<1,10
Z-POS-bdf800	1	800	<1,00
Z-POS-bdf1.5k	1	1.500	<1,10
Z-POS-bdf2k	1	2.000	<1,10
Z-POS-bdf2.5k	1	2.500	<1,10
Z-POS-bdf3k	1	3.000	<1,10
Z-POS-bdf4k	1	4.000	<1,10
Z-POS-bdf5k	1	5.000	<1,10
Z-POS-bdf7k	1	7.000	<1,10
Z-POS-bdf11k	1	11.000	<1,10
Z-POS-bdf13k	1	13.000	<1,10
Z-POS-bdf18k	1	18.000	<1,10
Z-POS-bdf20k	1	20.000	<1,10
Z-POS-bdf30k	1	30.000	<1,10
Z-POS-bdf35k	1	35.000	<1,10
Z-POS-bdf50k	1	50.000	<1,10
Z-POS-bdf60k	1	60.000	<1,00
Z-POS-bdf90k	1	90.000	<1,10
Z-POS-bdf150k	1	150.000	<1,10
Z-POS-bdf170k	1	170.000	<1,10
Z-POS-bdf200k	1	200.000	<1,10
Z-POS-bdf250k	1	250.000	<1,10
Z-POS-bdf300k	1	300.000	<1,10
Z-POS-bdf800k	1	800.000	<1,10
Z-POS-bdf1.5m	1	1.500.000	<1,10

Poly(butadiene d6)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-bdd200k	1	191.000	1,03

Poly(butylene terephthalate) broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pbt30k	1	30.000	<2,30
Z-POS-pbt40k	1	40.000	<2,30
Z-POS-pbt45k	1	45.000	<2,30
Z-POS-pbt50k	1	50.000	<2,30
Z-POS-pbt55k	1	55.000	<2,30

Poly(carbonate)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pc300	1	300	1
Z-POS-pc850	1	850	< 2,50
Z-POS-pc1k	1	1.000	< 2,50
Z-POS-pc1.5k	1	1.500	< 2,50
Z-POS-pc3k	1	3.000	< 2,50
Z-POS-pc4k	1	4.000	< 2,50
Z-POS-pc6k	1	6.000	< 2,50
Z-POS-pc7k	1	7.000	< 2,50
Z-POS-pc9k	1	9.000	< 2,50
Z-POS-pc10k	1	10.000	< 2,50
Z-POS-pc20k	1	20.000	< 2,50
Z-POS-pc25k	1	25.000	< 2,50
Z-POS-pc30k	1	30.000	< 2,50
Z-POS-pc50k	1	50.000	< 2,50
Z-POS-pc60k	1	60.000	< 2,90
Z-POS-pc90k	1	90.000	< 2,50
Z-POS-pc100k	1	100.000	< 2,90
Z-POS-pc130k	1	130.000	< 2,50
Z-POS-pc160k	1	160.000	3,21
Z-POS-pc195k	1	195.000	< 2,50

Poly(carbonate)-d4 deuterated

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pcd15k	1	15.000	2,62
Z-POS-pcd50k	1	50.000	1,97
Z-POS-pcd60k	1	60.000	1,90

Poly(dimethyl siloxane)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pdm162	1	162	1,00
Z-POS-pdm237	1	237	1,00
Z-POS-pdm311	1	311	1,00
Z-POS-pdm900	1	900	<1,40
Z-POS-pdm1.5k	1	1.500	<1,40
Z-POS-pdm2.5k	1	2.500	<1,40
Z-POS-pdm3.5k	1	3.500	<1,25
Z-POS-pdm4.5k	1	4.500	<1,25
Z-POS-pdm6k	1	6.000	<1,15
Z-POS-pdm10k	1	10.000	<1,15
Z-POS-pdm15k	1	15.000	<1,15
Z-POS-pdm25k	1	25.000	<1,15
Z-POS-pdm33k	1	33.000	<1,15
Z-POS-pdm40k	1	40.000	1,44
Z-POS-pdm50k	1	50.000	1,20
Z-POS-pdm75k	1	75.000	1,47
Z-POS-pdm110k	1	110.000	<1,15
Z-POS-pdm130k	1	130.000	<1,15
Z-POS-pdm150k	1	150.000	<1,15
Z-POS-pdm160k	1	160.000	<1,15
Z-POS-pdm200k	1	200.000	<1,15
Z-POS-pdm300k	1	300.000	1,11
Z-POS-pdm400k	1	400.000	1,17

Poly(dimethyl siloxane)-d6 deuterated

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pdd5k	1	5.000	1,14
Z-POS-pdd25k	1	25.000	1,12
Z-POS-pdd80k	1	80.000	1,39

Poly(ethyl methacrylate)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-em1.8k	1	1.800	<1,10
Z-POS-em4.5k	1	4.500	<1,10
Z-POS-em7.5k	1	7.500	<1,10
Z-POS-em45k	1	45.000	<1,10
Z-POS-em140k	1	140.000	<1,10
Z-POS-em165k	1	165.000	1,13
Z-POS-em270k	1	270.000	<1,10
Z-POS-em300k	1	300.000	<1,10
Z-POS-em350k	1	350.000	<1,10
Z-POS-em500k	1	500.000	<1,10

Poly(ethylene terephthalate) broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pet220	1	220	<3,00
Z-POS-pet2k	1	2.000	<3,00
Z-POS-pet5k	1	5.000	<3,40
Z-POS-pet10k	1	10.000	<3,20
Z-POS-pet13k	1	13.000	<3,70
Z-POS-pet20k	1	20.000	<2,30
Z-POS-pet25k	1	25.000	<2,40
Z-POS-pet30k	1	30.000	<2,80
Z-POS-pet50k	1	50.000	<2,80

Poly(ethylene) broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-peb50k	1	50.000	2,84
Z-POS-peb105k	1	105.000	10,10

Poly(iso-butylene)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pib112	1	112	1,00
Z-POS-pib224	1	224	1,00
Z-POS-pib350	1	350	<1,25
Z-POS-pib400	1	400	<1,25
Z-POS-pib700	1	700	<1,25
Z-POS-pib1.5k	1	1.500	<1,25
Z-POS-pib3.5k	1	3.500	<1,25
Z-POS-pib4k	1	4.000	<1,25
Z-POS-pib8k	1	8.000	<1,25
Z-POS-pib10k	1	10.000	<1,25
Z-POS-pib12k	1	12.000	<1,25
Z-POS-pib15k	1	15.000	<1,25
Z-POS-pib25k	1	25.000	<1,25
Z-POS-pib34k	1	34.000	<1,25
Z-POS-pib43k	1	43.000	<1,25
Z-POS-pib50k	1	50.000	<1,25
Z-POS-pib75k	1	75.000	<1,25
Z-POS-pib90k	1	90.000	<1,25
Z-POS-pib130k	1	130.000	<1,25
Z-POS-pib160k	1	160.000	<1,25
Z-POS-pib300k	1	300.000	<1,25
Z-POS-pib600k	1	600.000	1,29
Z-POS-pib700k	1	700.000	<1,25

Poly(iso-propyl methacrylate)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-ipm14k	1	14.000	<1,10
Z-POS-ipm20k	1	20.000	<1,10
Z-POS-ipm35k	1	35.000	<1,10
Z-POS-ipm45k	1	45.000	<1,10
Z-POS-ipm65k	1	65.000	<1,10
Z-POS-ipm160k	1	160.000	<1,10
Z-POS-ipm220k	1	220.000	<1,10
Z-POS-ipm400k	1	400.000	<1,10

Poly(isoprene 1,4)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pio600	1	600	<1,10
Z-POS-pio800	1	800	<1,10
Z-POS-pio900	1	900	1,17
Z-POS-pio1k	1	1.000	<1,10
Z-POS-pio1.3k	1	1.300	<1,10
Z-POS-pio1.5k	1	1.500	<1,10
Z-POS-pio2k	1	2.000	<1,10
Z-POS-pio4k	1	4.000	<1,10
Z-POS-pio4.5k	1	4.500	<1,10
Z-POS-pio7k	1	7.000	<1,10
Z-POS-pio8k	1	8.000	<1,10
Z-POS-pio10k	1	10.000	<1,10
Z-POS-pio13k	1	13.000	<1,10
Z-POS-pio15k	1	15.000	<1,10
Z-POS-pio17k	1	17.000	<1,10
Z-POS-pio20k	1	20.000	<1,10
Z-POS-pio23k	1	23.000	<1,10
Z-POS-pio35k	1	35.000	<1,10
Z-POS-pio50k	1	50.000	<1,10
Z-POS-pio60k	1	60.000	<1,10
Z-POS-pio70k	1	70.000	<1,10
Z-POS-pio80k	1	80.000	<1,10
Z-POS-pio100k	1	100.000	<1,10
Z-POS-pio150k	1	150.000	<1,10
Z-POS-pio200k	1	200.000	<1,10
Z-POS-pio270k	1	270.000	<1,10
Z-POS-pio300k	1	300.000	<1,10
Z-POS-pio500k	1	500.000	<1,10
Z-POS-pio600k	1	600.000	<1,10
Z-POS-pio700k	1	700.000	<1,10
Z-POS-pio1m	1	1.000.000	1,12
Z-POS-pio3m	1	3.000.000	1,23

Poly(isoprene 3,4)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pit900	1	900	<1,10
Z-POS-pit7k	1	7.000	<1,10
Z-POS-pit50k	1	50.000	<1,10
Z-POS-pit70k	1	70.000	1,12
Z-POS-pit250k	1	250.000	<1,10

Poly(methyl methacrylate)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-mm102	1	102	1,00
Z-POS-mm202	1	202	1,00
Z-POS-mm270	1	270	1,00
Z-POS-mm320	1	320	1,00
Z-POS-mm402	1	402	1,00
Z-POS-mm450	1	450	1,25
Z-POS-mm650	1	650	1,28
Z-POS-mm700	1	700	<1,10
Z-POS-mm850	1	850	1,44
Z-POS-mm2k	1	2.000	<1,10
Z-POS-mm2.5k	1	2.500	<1,10
Z-POS-mm3.5k	1	3.500	<1,10
Z-POS-mm4k	1	4.000	<1,10
Z-POS-mm5k	1	5.000	<1,10
Z-POS-mm8k	1	8.000	<1,10
Z-POS-mm10k	1	10.000	<1,10
Z-POS-mm13k	1	13.000	<1,10
Z-POS-mm15k	1	15.000	<1,10
Z-POS-mm20k	1	20.000	<1,10
Z-POS-mm23k	1	23.000	<1,10
Z-POS-mm30k	1	30.000	<1,10
Z-POS-mm40k	1	40.000	<1,10
Z-POS-mm45k	1	45.000	<1,10
Z-POS-mm60k	1	60.000	<1,10
Z-POS-mm65k	1	65.000	<1,10
Z-POS-mm80k	1	80.000	<1,10
Z-POS-mm90k	1	90.000	<1,10
Z-POS-mm100k	1	100.000	<1,10
Z-POS-mm130k	1	130.000	<1,10
Z-POS-mm150k	1	150.000	<1,10
Z-POS-mm160k	1	160.000	<1,10
Z-POS-mm170k	1	170.000	<1,10
Z-POS-mm180k	1	180.000	<1,10
Z-POS-mm200k	1	200.000	<1,10
Z-POS-mm250k	1	250.000	<1,10
Z-POS-mm280k	1	280.000	<1,10
Z-POS-mm350k	1	350.000	<1,10
Z-POS-mm500k	1	500.000	<1,10
Z-POS-mm600k	1	600.000	<1,10
Z-POS-mm650k	1	650.000	<1,10
Z-POS-mm700k	1	700.000	<1,10
Z-POS-mm750k	1	750.000	1,14
Z-POS-mm850k	1	850.000	<1,10
Z-POS-mm1.3m	1	1.300.000	1,11
Z-POS-mm1.5m	1	1.500.000	1,15
Z-POS-mm2.5M	1	2.500.000	1,16

Poly(methyl methacrylate) broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-mmb20k	1	20.000	2,80
Z-POS-mmb90k	1	90.000	1,80
Z-POS-mmb250k	1	250.000	2,20
Z-POS-mmb1.7m	1	1.700.000	2,22
Z-POS-mmb3m	1	3.000.000	1,96

Poly(methyl methacrylate) isotactic

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-mmi5k	1	5.000	<1,30
Z-POS-mmi12k	1	12.000	<1,30
Z-POS-mmi60k	1	60.000	<1,30

Poly(methyl methacrylate) syndiotactic

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-mms8k	1	8.000	<1,20
Z-POS-mms20k	1	20.000	1,30
Z-POS-mms25k	1	25.000	<1,20
Z-POS-mms50k	1	50.000	<1,20
Z-POS-mms400k	1	400.000	<1,20

Poly(methyl methacrylate)-d8 deuterated

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-mmd4k	1	4.000	1,12
Z-POS-mmd20k	1	20.000	<1,10
Z-POS-mmd35k	1	35.000	<1,10
Z-POS-mmd170k	1	170.000	1,31
Z-POS-mmd300k	1	315.000	1,64

Poly(n-butyl methacrylate)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-nb1k	1	1.000	1,11
Z-POS-nb1.5k	1	1.500	<1,10
Z-POS-nb3k	1	3.000	<1,10
Z-POS-nb4k	1	4.000	<1,10
Z-POS-nb5.5k	1	5.500	1,11
Z-POS-nb6.5k	1	6.500	1,11
Z-POS-nb8k	1	8.000	<1,10
Z-POS-nb12k	1	12.000	<1,10
Z-POS-nb20k	1	20.000	<1,10
Z-POS-nb45k	1	45.000	<1,10
Z-POS-nb60k	1	60.000	<1,10
Z-POS-nb80k	1	80.000	<1,10
Z-POS-nb110k	1	110.000	<1,10
Z-POS-nb220k	1	220.000	<1,10
Z-POS-nb330k	1	330.000	<1,10
Z-POS-nb500k	1	500.000	<1,10
Z-POS-nb600k	1	600.000	<1,10
Z-POS-nb700k	1	700.000	<1,10

Poly(n-propyl methacrylate)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-npm15k	1	15.000	<1,10
Z-POS-npm40k	1	40.000	1,50
Z-POS-npm50k	1	50.000	1,59
Z-POS-npm60k	1	60.000	<1,10

Poly(o-chloro styrene)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-ocl30k	1	30.000	1,90
Z-POS-ocl50k	1	50.000	<1,20
Z-POS-ocl100k	1	100.000	<1,20
Z-POS-ocl170k	1	170.000	<1,20
Z-POS-ocl300k	1	300.000	<1,20
Z-POS-ocl1.3m	1	1.300.000	1,76

Poly(p-methyl styrene)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-ppm60k	1	60.000	<1,10
Z-POS-ppm70k	1	70.000	<1,10
Z-POS-ppm130k	1	130.000	<1,10
Z-POS-ppm160k	1	160.000	<1,10

Poly(p-methyl styrene) deuterated

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pmd110k	1	110.000	<1,10

Poly(propylene) broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-ppb800	1	800	<4,00
Z-POS-ppb1.2k	1	1.200	<4,00
Z-POS-ppb6k	1	6.000	<4,00
Z-POS-ppb9k	1	9.000	<4,00
Z-POS-ppb45k	1	45.000	<4,00
Z-POS-ppb50k	1	50.000	<4,00
Z-POS-ppb60k	1	60.000	<4,00
Z-POS-ppb95k	1	95.000	<4,00
Z-POS-ppb105k	1	105.000	<4,00
Z-POS-ppb110k	1	110.000	<4,00
Z-POS-ppb135k	1	135.000	<4,00
Z-POS-ppb150k	1	150.000	<4,00
Z-POS-ppb180k	1	180.000	4,43
Z-POS-ppb200k	1	200.000	4,85
Z-POS-ppb210k	1	210.000	5,10
Z-POS-ppb230k	1	230.000	5,53
Z-POS-ppb270k	1	270.000	6,17
Z-POS-ppb300k	1	300.000	6,25
Z-POS-ppb350k	1	350.000	8,01

Poly(styrene)

Article No.	Amount [g]	Mw [D] \pm 10%	PDI
Z-POS-ps162	1	162	1,00
Z-POS-ps266	1	266	1,00
Z-POS-ps300	1	300	<1,10
Z-POS-ps370	1	370	1,00
Z-POS-ps500	1	500	<1,10
Z-POS-ps700	1	700	<1,10
Z-POS-ps800	1	800	<1,10
Z-POS-ps1k	1	1.000	<1,10
Z-POS-ps1.2k	1	1.200	<1,10
Z-POS-ps1.5k	1	1.500	<1,10
Z-POS-ps2k	1	2.000	<1,10
Z-POS-ps2.5k	1	2.500	<1,10
Z-POS-ps3.5k	1	3.500	<1,10
Z-POS-ps4k	1	4.000	<1,10
Z-POS-ps4.5k	1	4.500	<1,10
Z-POS-ps5k	1	5.000	<1,10
Z-POS-ps5.5k	1	5.500	<1,10
Z-POS-ps7k	1	7.000	<1,10
Z-POS-ps8k	1	8.000	<1,10
Z-POS-ps10k	1	10.000	<1,10
Z-POS-ps12k	1	12.000	<1,10
Z-POS-ps15k	1	15.000	<1,10
Z-POS-ps18k	1	18.000	<1,10
Z-POS-ps20k	1	20.000	<1,10
Z-POS-ps25k	1	25.000	<1,10
Z-POS-ps30k	1	30.000	<1,10
Z-POS-ps35k	1	35.000	<1,10
Z-POS-ps40k	1	40.000	<1,10
Z-POS-ps50k	1	50.000	<1,10
Z-POS-ps60k	1	60.000	<1,10
Z-POS-ps75k	1	75.000	<1,10
Z-POS-ps90k	1	90.000	<1,10
Z-POS-ps100k	1	100.000	<1,10
Z-POS-ps120k	1	120.000	<1,10
Z-POS-ps135k	1	135.000	<1,10
Z-POS-ps170k	1	170.000	<1,10
Z-POS-ps200k	1	200.000	<1,10
Z-POS-ps250k	1	250.000	<1,10
Z-POS-ps300k	1	300.000	<1,10
Z-POS-ps350k	1	350.000	<1,10
Z-POS-ps400k	1	400.000	<1,10
Z-POS-ps500k	1	500.000	<1,10
Z-POS-ps600k	1	600.000	<1,10
Z-POS-ps700k	1	700.000	<1,10
Z-POS-ps800k	1	800.000	<1,10
Z-POS-ps900k	1	900.000	<1,10
Z-POS-ps1m	1	1.000.000	<1,10
Z-POS-ps1.1m	1	1.100.000	1,14
Z-POS-ps1.3m	1	1.300.000	<1,10
Z-POS-ps1.5m	1	1.500.000	<1,10
Z-POS-ps1.8m	1	1.800.000	<1,10
Z-POS-ps2.5m	1	2.500.000	<1,10
Z-POS-ps3m	1	3.000.000	<1,30

Poly(styrene)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-ps6m	1	6.000.000	<1,30
Z-POS-ps8m	1	8.000.000	<1,30
Z-POS-ps11m	1	11.000.000	<1,30
Z-POS-ps23m	1	23.000.000	<1,30

Poly(styrene) DIN cert.

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-dps162	1	162	1,00
Z-POS-dps300	1	300	<1,20
Z-POS-dps700	1	700	<1,20
Z-POS-dps1.1k	1	1.100	<1,20
Z-POS-dps1.5k	1	1.500	<1,20
Z-POS-dps2.5k	1	2.500	<1,05
Z-POS-dps3k	1	3.000	<1,05
Z-POS-dps4k	1	4.000	<1,05
Z-POS-dps5k	1	5.000	<1,05
Z-POS-dps6k	1	6.000	<1,05
Z-POS-dps8k	1	8.000	<1,05
Z-POS-dps20k	1	20.000	<1,05
Z-POS-dps30k	1	30.000	<1,05
Z-POS-dps50k	1	50.000	<1,05
Z-POS-dps80k	1	80.000	<1,05
Z-POS-dps100k	1	100.000	<1,05
Z-POS-dps120k	1	120.000	<1,05
Z-POS-dps130k	1	130.000	<1,05
Z-POS-dps200k	1	200.000	<1,05
Z-POS-dps250k	1	250.000	<1,05
Z-POS-dps600k	1	600.000	<1,05
Z-POS-dps800k	1	800.000	<1,05
Z-POS-dps1m	1	1.000.000	<1,20
Z-POS-dps2m	1	2.000.000	<1,20
Z-POS-dps4m	1	4.000.000	<1,20
Z-POS-dps5m	1	5.000.000	<1,20
Z-POS-dps7m	1	7.000.000	<1,20

Poly(styrene) broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-psb25k	1	25.000	<3,00
Z-POS-psb50k	1	50.000	<3,00
Z-POS-psb90k	1	90.000	<3,00
Z-POS-psb130k	1	130.000	<3,00
Z-POS-psb170k	1	170.000	<3,00
Z-POS-psb200k	1	200.000	<3,00
Z-POS-psb250k	1	250.000	<3,00
Z-POS-psb330k	1	330.000	<3,00
Z-POS-psb350k	1	350.000	<3,00
Z-POS-psb400k	1	400.000	<3,00
Z-POS-psb500k	1	500.000	<3,00
Z-POS-psb600k	1	600.000	<3,00
Z-POS-psb700k	1	700.000	<3,00
Z-POS-psb3m	1	3.000.000	<3,00

Poly(styrene)-d8 deuterated

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-psd900	1	900	<1,10
Z-POS-psd3k	1	3.000	<1,10
Z-POS-psd4k	1	4.000	<1,10
Z-POS-psd5k	1	5.000	<1,10
Z-POS-psd6k	1	6.000	<1,10
Z-POS-psd7k	1	7.000	<1,10
Z-POS-psd12k	1	12.000	<1,10
Z-POS-psd15k	1	15.000	<1,10
Z-POS-psd25k	1	25.000	<1,10
Z-POS-psd30k	1	30.000	<1,10
Z-POS-psd40k	1	40.000	<1,10
Z-POS-psd50k	1	50.000	<1,10
Z-POS-psd70k	1	70.000	<1,10
Z-POS-psd80k	1	80.000	<1,10
Z-POS-psd90k	1	90.000	<1,10
Z-POS-psd100k	1	100.000	<1,10
Z-POS-psd110k	1	110.000	<1,10
Z-POS-psd120k	1	120.000	<1,10
Z-POS-psd130k	1	130.000	<1,10
Z-POS-psd160k	1	160.000	<1,10
Z-POS-psd170k	1	170.000	<1,10
Z-POS-psd250k	1	250.000	<1,10
Z-POS-psd300k	1	300.000	<1,10
Z-POS-psd400k	1	400.000	<1,10
Z-POS-psd500k	1	500.000	<1,10
Z-POS-psd600k	1	600.000	<1,10
Z-POS-psd800k	1	800.000	<1,10
Z-POS-psd1m	1	1.000.000	<1,10
Z-POS-psd1.2m	1	1.200.000	<1,10
Z-POS-psd1.5m	1	1.500.000	<1,10
Z-POS-psd1.7m	1	1.700.000	<1,10

Poly(t-butyl acrylate)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-tba30k	1	30.000	<1,25
Z-POS-tba40k	1	40.000	<1,25
Z-POS-tba50k	1	50.000	<1,25
Z-POS-tba55k	1	55.000	<1,25
Z-POS-tba60k	1	60.000	<1,25
Z-POS-tba70k	1	70.000	<1,25
Z-POS-tba90k	1	90.000	<1,25
Z-POS-tba100k	1	100.000	<1,10
Z-POS-tba150k	1	150.000	<1,10
Z-POS-tba180k	1	180.000	<1,10
Z-POS-tba220k	1	220.000	<1,10
Z-POS-tba285k	1	285.000	<1,10
Z-POS-tba390k	1	390.000	<1,10
Z-POS-tba440k	1	440.000	<1,10
Z-POS-tba600k	1	600.000	<1,10
Z-POS-tba900k	1	900.000	<1,10
Z-POS-tba1.2m	1	1.200.000	<1,25
Z-POS-tba3m	1	3.000.000	<1,25

Poly(tert.butylmethacrylate)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-tb600	1	600	<1,25
Z-POS-tb1.5k	1	1.500	<1,25
Z-POS-tb1.7k	1	1.700	<1,25
Z-POS-tb2k	1	2.000	<1,25
Z-POS-tb2.3k	1	2.300	<1,25
Z-POS-tb4k	1	4.000	<1,25
Z-POS-tb5k	1	5.000	<1,25
Z-POS-tb6.5k	1	6.500	<1,25
Z-POS-tb8k	1	8.000	<1,25
Z-POS-tb10k	1	10.000	<1,10
Z-POS-tb15k	1	15.000	<1,10
Z-POS-tb20k	1	20.000	<1,10
Z-POS-tb25k	1	25.000	<1,10
Z-POS-tb30k	1	30.000	<1,10
Z-POS-tb35k	1	35.000	<1,10
Z-POS-tb40k	1	40.000	<1,10
Z-POS-tb45k	1	45.000	<1,10
Z-POS-tb60k	1	60.000	<1,10
Z-POS-tb100k	1	100.000	<1,10
Z-POS-tb135k	1	135.000	<1,10
Z-POS-tb200k	1	200.000	<1,10
Z-POS-tb250k	1	250.000	<1,10
Z-POS-tb280k	1	280.000	<1,10
Z-POS-tb310k	1	310.000	<1,10
Z-POS-tb550k	1	550.000	<1,15
Z-POS-tb650k	1	650.000	<1,15
Z-POS-tb850k	1	850.000	<1,15
Z-POS-tb1.2m	1	1.200.000	1,47

Poly(tetrahydrofurane)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-ptb90	1	90	1,00
Z-POS-ptb230	1	230	<1,60
Z-POS-ptb300	1	300	<1,60
Z-POS-ptb900	1	900	<1,60
Z-POS-ptb1.5k	1	1.500	2,01
Z-POS-ptb5k	1	5.000	2,19
Z-POS-ptb9k	1	9.000	<1,60
Z-POS-ptb12k	1	12.000	<1,60
Z-POS-ptb15k	1	15.000	<1,60
Z-POS-ptb20k	1	20.000	<1,60
Z-POS-ptb50k	1	50.000	<1,60
Z-POS-ptb90k	1	90.000	<1,60
Z-POS-ptb110k	1	110.000	<1,60

Poly(vinyl acetate) broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pva17k	1	17.000	<3,00
Z-POS-pva50k	1	50.000	<3,00
Z-POS-pva120k	1	120.000	<3,00
Z-POS-pva150k	1	150.000	<3,00
Z-POS-pva190k	1	190.000	<3,00
Z-POS-pva200k	1	200.000	<3,00
Z-POS-pva250k	1	250.000	<3,00
Z-POS-pva300k	1	300.000	4,19

Poly(vinyl butyral) broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pvb70k	1	70.000	<4,00
Z-POS-pvb90k	1	90.000	<4,00
Z-POS-pvb130k	1	130.000	<4,00
Z-POS-pvb160k	1	160.000	4,03
Z-POS-pvb170k	1	170.000	<4,00

Poly(vinyl chloride)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pvc127	1	127	1,00
Z-POS-pvc40k	1	40.000	<1,20
Z-POS-pvc50k	1	50.000	<1,20
Z-POS-pvc60k	1	60.000	<1,20
Z-POS-pvc70k	1	70.000	<1,20
Z-POS-pvc100k	1	100.000	<1,35
Z-POS-pvc120k	1	120.000	<1,35
Z-POS-pvc150k	1	150.000	<1,35
Z-POS-pvc170k	1	170.000	<1,35
Z-POS-pvc220k	1	220.000	<1,35
Z-POS-pvc300k	1	300.000	1,40

Polyethylene

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pe1k	1	1.000	1.11
Z-POS-pe2k	1	2.000	1.11
Z-POS-pe16k	1	16.000	<1.60
Z-POS-pe20k	1	20.000	<1.60
Z-POS-pe25k	1	25.000	<1.60
Z-POS-pe30k	1	30.000	1.84
Z-POS-pe35k	1	35.000	<1.60
Z-POS-pe60k	1	60.000	<1.60
Z-POS-pe75k	1	75.000	<1.60
Z-POS-pe85k	1	85.000	<1.60
Z-POS-pe90k	1	90.000	<1.60
Z-POS-pe115k	1	115.000	<1.60
Z-POS-pe150k	1	150.000	<1.60
Z-POS-pe170k	1	170.000	<1.60
Z-POS-pe180k	1	180.000	<1.60

Polypropylene

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pp700	0.25	700	1.06
Z-POS-pp800	0.25	800	1.05
Z-POS-pp900	0.25	900	1.04
Z-POS-pp1k	0.25	1.000	1.05
Z-POS-pp1.1k	0.25	1.100	1.03
Z-POS-pp1.2k	0.25	1.200	1.05
Z-POS-pp1.3k	0.25	1.300	1.05
Z-POS-pp1.4k	0.25	1.400	1.03
Z-POS-pp1.5k	0.25	1.500	1.05
Z-POS-pp1.7k	0.25	1.700	1.05
Z-POS-pp1.8k	0.25	1.800	1.09
Z-POS-pp1.9k	0.25	1.900	1.06
Z-POS-pp2.1k	0.25	2.100	1.06
Z-POS-pp2.3k	0.25	2.300	1.03
Z-POS-pp2.6k	0.25	2.600	1.05
Z-POS-pp2.8k	0.25	2.800	1.06

Aqueous Polymers

Dextrane

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-dxt180	1	180	1,00
Z-POS-dxt342	1	342	1,00
Z-POS-dxt504	1	504	1,00
Z-POS-dxt1.3k	1	1.300	<1,20
Z-POS-dxt1.5k	1	1.500	<1,20
Z-POS-dxt5k	1	5.000	<1,80
Z-POS-dxt12k	1	12.000	<1,80
Z-POS-dxt25k	1	25.000	<1,80
Z-POS-dxt50k	1	50.000	<1,80
Z-POS-dxt80k	1	80.000	<1,80
Z-POS-dxt150k	1	150.000	<1,80
Z-POS-dxt270k	1	270.000	<1,80
Z-POS-dxt410k	1	410.000	<1,80
Z-POS-dxt670k	1	670.000	<1,80
Z-POS-dxt3.8m	1	3.800.000	2,53

Dextrane DIN cert.

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-ddx180	1	180	1,00
Z-POS-ddx342	1	342	1,00
Z-POS-ddx504	1	504	1,00
Z-POS-ddx1.3k	1	1.300	<1,20
Z-POS-ddx1.5k	1	1.500	<1,20
Z-POS-ddx5k	1	5.000	<1,80
Z-POS-ddx12k	1	12.000	<1,80
Z-POS-ddx25k	1	25.000	<1,80
Z-POS-ddx50k	1	50.000	<1,80
Z-POS-ddx80k	1	80.000	<1,80
Z-POS-ddx150k	1	150.000	<1,80
Z-POS-ddx270k	1	270.000	<1,80
Z-POS-ddx410k	1	410.000	<1,80
Z-POS-ddx670k	1	670.000	<1,80
Z-POS-ddx3.8m	1	3.800.000	2,53

Dextrane broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-dxb4k	1	4.000	<2,10
Z-POS-dxb10k	1	11.000	<2,10
Z-POS-dxb40k	1	40.000	<2,10
Z-POS-dxb70k	1	70.000	<2,10
Z-POS-dxb500k	1	500.000	<2,10
Z-POS-dxb2m	1	2.000.000	<2,10
Z-POS-dxb17m	1	17.000.000	<2,10

Hydroxy-ethyl starch

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-hes25k	1	25.000	<2,30
Z-POS-hes60k	1	60.000	<2,30
Z-POS-hes120k	1	120.000	<2,30
Z-POS-hes140k	1	140.000	2,58
Z-POS-hes160k	1	160.000	3,83
Z-POS-hes400k	1	400.000	9,81
Z-POS-hes1m	1	1.000.000	<2,30
Z-POS-hes1.8m	1	1.800.000	<2,30
Z-POS-hes2m	1	2.000.000	<2,30

Hydroxy-propyl cellulose

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-hpc12k	1	12.000	<2,00
Z-POS-hpc15k	1	15.000	<2,00
Z-POS-hpc20k	1	20.000	<2,00
Z-POS-hpc30k	1	30.000	<2,00
Z-POS-hpc33k	1	33.000	<2,00
Z-POS-hpc40k	1	40.000	<2,00
Z-POS-hpc60k	1	60.000	<2,00

Hydroxy-propyl cellulose broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-hpb30k	1	30.000	<6,00
Z-POS-hpb50k	1	50.000	<6,00
Z-POS-hpb55k	1	55.000	<6,00
Z-POS-hpb60k	1	60.000	<6,00
Z-POS-hpb70k	1	70.000	<6,00
Z-POS-hpb80k	1	80.000	<6,00
Z-POS-hpb90k	1	90.000	<6,00
Z-POS-hpb100k	1	100.000	<6,00
Z-POS-hpb150k	1	150.000	<6,00
Z-POS-hpb180k	1	180.000	<6,00
Z-POS-hpb200k	1	200.000	<6,00
Z-POS-hpb300k	1	300.000	7,44
Z-POS-hpb400k	1	400.000	<6,00
Z-POS-hpb500k	1	500.000	<6,00
Z-POS-hpb550k	1	550.000	<6,00
Z-POS-hpb600k	1	600.000	<6,00
Z-POS-hpb700k	1	700.000	<6,00
Z-POS-hpb900k	1	900.000	<6,00

Poly(2-vinyl pyridine) quaternated

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pvq4k	1	4.000	<1,20
Z-POS-pvq7k	1	7.000	<1,20
Z-POS-pvq18k	1	18.000	<1,20
Z-POS-pvq30k	1	30.000	<1,20
Z-POS-pvq50k	1	50.000	<1,20
Z-POS-pvq130k	1	130.000	<1,20
Z-POS-pvq300k	1	300.000	<1,20
Z-POS-pvq700k	1	700.000	<1,20

Poly(acrylamid) broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pam3k	1	3.000	<2,60
Z-POS-pam12k	1	12.000	<2,60
Z-POS-pam20k	1	20.000	<2,60
Z-POS-pam60k	1	60.000	<2,60
Z-POS-pam80k	1	80.000	<2,60
Z-POS-pam400k	1	400.000	<2,60
Z-POS-pam500k	1	500.000	<2,60
Z-POS-pam1.1m	1	1.100.000	<2,60
Z-POS-pam6m	1	6.000.000	<2,60
Z-POS-pam9m	1	9.000.000	<2,60

Poly(acrylic acid) sodium salt

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-paa2k	1	2.000	<1,80
Z-POS-paa2.5k	1	2.500	<1,80
Z-POS-paa4k	1	4.000	<1,80
Z-POS-paa6k	1	6.000	<1,80
Z-POS-paa8k	1	8.000	<1,80
Z-POS-paa20k	1	20.000	<1,80
Z-POS-paa40k	1	40.000	<1,80
Z-POS-paa80k	1	80.000	<1,80
Z-POS-paa130k	1	130.000	<1,80
Z-POS-paa150k	1	150.000	<1,80
Z-POS-paa250k	1	250.000	<1,80
Z-POS-paa350k	1	350.000	<1,80
Z-POS-paa500k	1	500.000	<1,80
Z-POS-paa600k	1	600.000	<1,80
Z-POS-paa700k	1	700.000	<1,80
Z-POS-paa800k	1	800.000	<1,80
Z-POS-paa1m	1	1.000.000	<1,80
Z-POS-paa1.4m	1	1.400.000	<1,80

Poly(DADMAC)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-dad45k	1	45.000	<2,00
Z-POS-dad150k	1	150.000	<2,00
Z-POS-dad200k	1	200.000	<2,00
Z-POS-dad300k	1	300.000	2,26
Z-POS-dad450k	1	450.000	2,77

Poly(allyl amine) hydrochloride

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pah35k	1	35.000	<1,75
Z-POS-pah100k	1	100.000	<1,75

Poly(ethylene glycol)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-peg62	1	62	1,00
Z-POS-peg106	1	106	1,00
Z-POS-peg150	1	150	1,00
Z-POS-peg194	1	194	1,00
Z-POS-peg238	1	238	1,00
Z-POS-peg282	1	282	1,00
Z-POS-peg330	1	330	<1,10
Z-POS-peg400	1	400	<1,10
Z-POS-peg600	1	600	<1,10
Z-POS-peg1k	1	1.000	<1,10
Z-POS-peg1.5k	1	1.500	<1,10
Z-POS-peg2k	1	2.000	<1,10
Z-POS-peg3k	1	3.000	<1,10
Z-POS-peg4k	1	4.000	<1,10
Z-POS-peg6k	1	6.000	<1,10
Z-POS-peg8k	1	8.000	<1,10
Z-POS-peg11k	1	11.000	<1,10
Z-POS-peg14k	1	14.000	<1,10
Z-POS-peg23k	1	23.000	<1,10
Z-POS-peg40k	1	40.000	1,14

Poly(ethylene glycol) DIN cert.

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-dpg62	1	62	1,00
Z-POS-dpg106	1	106	1,00
Z-POS-dpg150	1	150	1,00
Z-POS-dpg194	1	194	1,00
Z-POS-dpg238	1	238	1,00
Z-POS-dpg282	1	282	1,00
Z-POS-dpg330	1	330	<1,10
Z-POS-dpg400	1	400	<1,10
Z-POS-dpg600	1	600	<1,10
Z-POS-dpg1k	1	1.000	<1,10
Z-POS-dpg1.5k	1	1.500	<1,10
Z-POS-dpg2k	1	2.000	<1,10
Z-POS-dpg3k	1	3.000	<1,10
Z-POS-dpg4k	1	4.000	<1,10
Z-POS-dpg6k	1	6.000	<1,10
Z-POS-dpg8k	1	8.000	<1,10
Z-POS-dpg11k	1	11.000	<1,10
Z-POS-dpg14k	1	14.000	<1,10
Z-POS-dpg23k	1	23.000	<1,10
Z-POS-dpg40k	1	40.000	1,14

Poly(ethylene oxide)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-peo20k	1	20.000	<1,25
Z-POS-peo55k	1	55.000	<1,25
Z-POS-peo70k	1	70.000	<1,25
Z-POS-peo90k	1	90.000	<1,25
Z-POS-peo120k	1	120.000	<1,25
Z-POS-peo130k	1	130.000	1,33
Z-POS-peo180k	1	180.000	<1,25
Z-POS-peo200k	1	200.000	<1,25
Z-POS-peo220k	1	220.000	<1,25
Z-POS-peo240k	1	240.000	<1,25
Z-POS-peo300k	1	300.000	1,26
Z-POS-peo400k	1	400.000	<1,25
Z-POS-peo450k	1	450.000	<1,25
Z-POS-peo500k	1	500.000	<1,25
Z-POS-peo600k	1	600.000	<1,25
Z-POS-peo1.5m	1	1.500.000	<1,25

Poly(ethylene oxide)-d4 deuterated

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-ped8k	1	8.000	1,06
Z-POS-ped13k	1	13.000	1,08
Z-POS-ped500k	1	460.000	1,44

Poly(methacrylic acid) sodium salt

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pma1.3k	1	1.300	<1,20
Z-POS-pma4k	1	4.000	<1,20
Z-POS-pma5k	1	5.000	<1,20
Z-POS-pma6k	1	6.000	<1,20
Z-POS-pma8k	1	8.000	<1,10
Z-POS-pma15k	1	15.000	<1,10
Z-POS-pma 20k	1	20.000	<1,10
Z-POS-pma30k	1	30.000	<1,10
Z-POS-pma40k	1	40.000	<1,10
Z-POS-pma60k	1	60.000	<1,10
Z-POS-pma75k	1	75.000	<1,10
Z-POS-pma130k	1	130.000	<1,10
Z-POS-pma200k	1	200.000	<1,10
Z-POS-pma350k	1	350.000	<1,15
Z-POS-pma400k	1	400.000	<1,15
Z-POS-pma500k	1	500.000	<1,15
Z-POS-pma800k	1	800.000	<1,15

Poly(methacrylic acid)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pms 6k	1	6.000	<1,10
Z-POS-pms12k	1	12.000	<1,10
Z-POS-pms15k	1	15.000	<1,10
Z-POS-pms60k	1	60.000	<1,10
Z-POS-pms80k	1	80.000	<1,10
Z-POS-pms110k	1	110.000	<1,10
Z-POS-pms200k	1	200.000	<1,10
Z-POS-pms350k	1	350.000	1,11
Z-POS-pms500	1	500.000	<1,10
Z-POS-pms800k	1	800.000	1,47

Poly(propylene glycol)

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-ppg76	1	76	1,00
Z-POS-ppg134	1	134	1,00
Z-POS-ppg192	1	192	1,00
Z-POS-ppg450	1	450	<1,10
Z-POS-ppg800	1	800	<1,10
Z-POS-ppg1k	1	1.000	<1,10
Z-POS-ppg2k	1	2.000	<1,10
Z-POS-ppg3k	1	3.000	<1,10
Z-POS-ppg5k	1	5.000	<1,10

Poly(styrene sulfonate) sodium salt

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pss208	1	208	1,00
Z-POS-pss800	1	800	<1,20
Z-POS-pss1.3k	1	1.300	<1,20
Z-POS-pss2.2k	1	2.200	<1,20
Z-POS-pss4.5k	1	4.500	<1,20
Z-POS-pss5k	1	5.000	<1,20
Z-POS-pss6.5k	1	6.500	<1,20
Z-POS-pss10k	1	10.000	<1,20
Z-POS-pss13k	1	13.000	<1,20
Z-POS-pss15k	1	15.000	<1,20
Z-POS-pss20k	1	20.000	<1,20
Z-POS-pss30k	1	30.000	<1,20
Z-POS-pss45k	1	45.000	<1,20
Z-POS-pss80k	1	80.000	<1,20
Z-POS-pss140k	1	140.000	<1,20
Z-POS-pss150k	1	150.000	<1,20
Z-POS-pss200k	1	200.000	<1,20
Z-POS-pss350k	1	350.000	<1,20
Z-POS-pss900k	1	900.000	<1,20
Z-POS-pss1.3m	1	1.300.000	<1,20
Z-POS-pss2.5m	1	2.500.000	<1,20
Z-POS-pss3.5m	1	3.500.000	<1,20
Z-POS-pss6m	1	6.000.000	<1,20

Poly(vinyl alcohol) broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pvo6k	1	6.000	<3,00
Z-POS-pvo12k	1	12.000	<3,00
Z-POS-pvo30k	1	30.000	<3,00
Z-POS-pvo35k	1	35.000	<3,00
Z-POS-pvo40k	1	40.000	<3,00
Z-POS-pvo70k	1	70.000	<3,00
Z-POS-pvo80k	1	80.000	<3,00
Z-POS-pvo100k	1	100.000	<3,00
Z-POS-pvo110k	1	110.000	<3,00
Z-POS-pvo160k	1	160.000	<3,00

Poly(vinyl pyrrolidone) broad standard

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-vpr15k	1	15.000	10,87
Z-POS-vpr17k	1	17.000	<9,00
Z-POS-vpr20k	1	20.000	<9,00
Z-POS-vpr70k	1	70.000	<9,00
Z-POS-vpr80k	1	80.000	<9,00
Z-POS-vpr100k	1	100.000	<9,00
Z-POS-vpr360k	1	360.000	<9,00
Z-POS-vpr1m	1	1.000.000	<9,00
Z-POS-vpr1.3m	1	1.300.000	<9,00
Z-POS-vpr1.5m	1	1.500.000	<9,00
Z-POS-vpr3.5m	1	3.500.000	10,57

Pullulane

Article No.	Amount [g]	Mw [D]±10%	PDI
Z-POS-pul180	1	180	1,00
Z-POS-pul342	1	342	1,00
Z-POS-pul504	1	504	1,00
Z-POS-pul1.5k	1	1.500	1,51