

# PN7140

FFF Solvent Organizer Module



# PN7140 Solvent Organizer

## Features

### Operation Principle:

Holder for bottles and possibility of direct connection of up to 8 eluent lines via the integrated PEEK connector block.

### Bottle Volumes:

For two 2 L, four 1 L, six 0.5 L or eight 0.25 L eluent/waste bottles.

Special versions may be provided by Postnova as well. Please inquire for possible options, which will depend on the requested customer application.

The PN7140 will be delivered including 2 x 1 liter bottles with bottle caps and 1/8" OD tubing.

## Ordering Information

S-REA-7140-001      2-Channel Version

Extendable up to 8 channels



## Specifications

- Working Principle:  
By the use of a special collection tray
- Flow Rate:  
Max. 10 mL/min for each channel
- Channels:  
Standard: 2 channels  
Extendable up to 8 channels
- Wetted Parts:  
PEEK, ETFE, PFA, UHMW-PE
- Dimensions:  
Width x Height x Length  
270 x 90 x 430 mm
- Weight:  
4 kg

## Contact

- Postnova Analytics GmbH  
86899 Landsberg, GERMANY  
T: +49 8191 985 688 0
- Postnova Analytics UK Ltd.  
Worcestershire, WR14 3SZ, UK  
T: +44 1684 585167
- Postnova Analytics Inc.  
Salt Lake City, UT 84102, USA  
T: +1 801 521 2004
- Postnova North Europe  
01630 Vantaa, FINLAND  
T: +358 9 8545 510

[info@postnova.com](mailto:info@postnova.com)  
[www.postnova.com](http://www.postnova.com)