

## Membranes for Flow Field-Flow Fractionation

Postnova offers a broad range of special membranes for Flow FFF Systems. Please check the system you use and select the corresponding membrane from the list below or contact our customer support.

<b>AF2000 MF Series (Asymmetric Flow FFF) Preparative Version</b>		<b>Qty / Pack</b>
<b>AF4-Membranes for Aqueous/Organic Solvents</b>		
Z-AF4-MEM-601-300D	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 300 Da cut-off; Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-601-1KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 1 kDa cut-off; Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-601-4KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 4 kDa cut-off; Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-601-5KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 5 kDa cut-off; Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-601-10KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off; Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-601-20KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off; Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-601-50KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 50 kDa cut-off; Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-601-150KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off; Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-602-1KD	NovaRC AQU Membrane pre-cut for 700 Channel Series, aqueous, 1 kDa cut-off; reg. Cellulose (RC), Size 110 x 294 mm	1
Z-AF4-MEM-602-5KD	NovaRC AQU/ORG Membrane pre-cut for 700 Channel Series, aqueous & organic, 5 kDa cut-off; reg. Cellulose (RC), Size 110 x 294 mm	1
Z-AF4-MEM-602-10KD	NovaRC AQU/ORG Membrane pre-cut for 700 Channel Series, aqueous & organic, 10 kDa cut-off; reg. Cellulose (RC), Size 110 x 294 mm	1
Z-AF4-MEM-602-30KD	NovaRC AQU Membrane pre-cut for 700 Channel Series, aqueous, 30 kDa cut-off; reg. Cellulose (RC), Size 110 x 294 mm	1
Z-AF4-MEM-602-100KD	NovaRC AQU Membrane pre-cut for 700 Channel Series, aqueous, 100 kDa cut-off; reg. Cellulose (RC), Size 110 x 294 mm	1
Z-AF4-MEM-605-10KD	NovaRC AMPH amphiphilic Membrane pre-cut for 700 Channel Series, aqueous, 10kDa cut-off, reg.Cellulose, for amphiphilic polymers/cationic polyelectrolytes, Size 110 x 294 mm	1
Z-AF4-MEM-603-10KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Cellulose Triacetate (CTA), Size 110 x 294 mm	1
Z-AF4-MEM-603-20KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off, Cellulose Triacetate (CTA), Size 110 x 294 mm	1
Z-AF4-MEM-604-150KD	NovaTF AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off, PVDF hydrophilic, Size 110 x 294 mm	1
Z-AF4-MEM-607-1KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, aqueous, 1 kDa cut-off, Polyacrylnitrile (PAN), Size 110 x 294 mm	1
Z-AF4-MEM-607-3KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, aqueous, 3 kDa cut-off, Polyacrylnitrile (PAN), Size 110 x 294 mm	1
Z-AF4-MEM-607-5KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, aqueous, 5 Da cut-off, Polyacrylnitrile (PAN), Size 110 x 294 mm	1
Z-AF4-MEM-607-10KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Polyacrylnitrile (PAN), Size 110 x 29 mm	1
Z-AF4-MEM-607-30KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, aqueous, 30 kDa cut-off, Polyacrylnitrile (PAN), Size 110 x 294 mm	1

<b>AF2000 MF Series (Asymmetric Flow FFF) Semipreparative Version</b>		<b>Qty / Pack</b>
<b>AF4-Membranes for Aqueous/Organic Solvents</b>		
Z-AF4-MEM-741-300D	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 300 Da cut-off; Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-741-1KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 1 kDa cut-off; Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-741-4KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 4 kDa cut-off; Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-741-5KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 5 kDa cut-off; Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-741-10KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-741-20KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off, Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-741-50KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 50 kDa cut-off, Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-741-150KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off, Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-742-1KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 1 kDa cut-off; reg. Cellulose (RC), Size 56 x 283 mm	1
Z-AF4-MEM-742-5KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 5 kDa cut-off; reg. Cellulose (RC), Size 56 x 283 mm	1
Z-AF4-MEM-742-10KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 10 kDa cut-off; reg. Cellulose (RC), Size 56 x 283 mm	1
Z-AF4-MEM-742-30KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 30 kDa cut-off; reg. Cellulose (RC), Size 56 x 283 mm	1
Z-AF4-MEM-742-100KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 100 kDa cut-off; reg. Cellulose (RC), Size 56 x 283 mm	1
Z-AF4-MEM-745-10KD	NovaRC AMPH amphiphilic Membrane pre-cut for 700 Channel Series, aqueous, 10kDa cut-off, reg.Cellulose, for amphiphilic polymers/cationic polyelectrolytes, Size 56 x 283 mm	1
Z-AF4-MEM-743-10KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Cellulose Triacetate (CTA), Size 56 x 283 mm	1
Z-AF4-MEM-743-20KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off, Cellulose Triacetate (CTA), Size 56 x 283 mm	1
Z-AF4-MEM-744-150KD	NovaTF AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off, PVDF hydrophilic, Size 56 x 283 mm	1
Z-AF4-MEM-747-1KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, aqueous, 1 kDa cut-off, Polyacrylnitrile (PAN), Size 56 x 283 mm	1
Z-AF4-MEM-747-3KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, aqueous, 3 kDa cut-off, Polyacrylnitrile (PAN), Size 56 x 283 mm	1
Z-AF4-MEM-747-5KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, aqueous, 5 Da cut-off, Polyacrylnitrile (PAN), Size 56 x 283 mm	1
Z-AF4-MEM-747-10KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Polyacrylnitrile (PAN), Size 56 x 283 mm	1
Z-AF4-MEM-747-30KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, aqueous, 30 kDa cut-off, Polyacrylnitrile (PAN), Size 56 x 283 mm	1

<b>AF2000 MF Series (Asymmetric Flow FFF) Analytical Version</b>		<b>Qty / Pack</b>
<b>AF4-Membranes for Aqueous/Organic Solvents</b>		
Z-AF4-MEM-711-300D	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 300 Da cut-off; Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-711-1KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 1 kDa cut-off; Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-711-4KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 4 kDa cut-off; Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-711-5KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 5 kDa cut-off; Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-711-10KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-711-20KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-711-50KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 50 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-711-150KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-712-1KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 1 kDa cut-off; reg. Cellulose (RC), Size 25 x 282 mm	10
Z-AF4-MEM-712-5KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 5 kDa cut-off; reg. Cellulose (RC), Size 25 x 282 mm	10
Z-AF4-MEM-712-10KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 10 kDa cut-off; reg. Cellulose (RC), Size 25 x 282 mm	10
Z-AF4-MEM-712-30KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 30 kDa cut-off; reg. Cellulose (RC), Size 25 x 282 mm	10
Z-AF4-MEM-712-100KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 100 kDa cut-off; reg. Cellulose (RC), Size 25 x 282 mm	10
Z-AF4-MEM-715-10KD	NovaRC AMPH amphiphilic Membrane pre-cut for 700 Channel Series, aqueous, 10kDa, reg.Cellulose, for amphiphilic polymers/cationic polyelectrolytes, Size 25 x 282 mm	10
Z-AF4-MEM-713-10KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Cellulose Triacetate (CTA), Size 25 x 282 mm	10
Z-AF4-MEM-713-20KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off, Cellulose Triacetate (CTA), Size 25 x 282 mm	10
Z-AF4-MEM-714-50KD	NovaTF AQU Membrane pre-cut for 700 Channel Series, aqueous, 50 kDa cut-off, PVDF hydrophilic, Size 25 x 282 mm	10
Z-AF4-MEM-714-150KD	NovaTF AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off, PVDF hydrophilic, Size 25 x 282 mm	10
Z-AF4-MEM-717-1KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, 1 kDa Polyacrylonitrile (PAN) Membrane; aqueous; Size: 25 x 282 mm	10
Z-AF4-MEM-717-3KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, 3 kDa Polyacrylonitrile (PAN) Membrane; aqueous; Size: 25 x 282 mm	10
Z-AF4-MEM-717-5KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, 5 kDa Polyacrylonitrile (PAN) Membrane; aqueous; Size: 25 x 282 mm	10
Z-AF4-MEM-717-10KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, 10 kDa Polyacrylonitrile (PAN) Membrane; aqueous; Size: 25 x 282 mm	10
Z-AF4-MEM-717-30KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, 30 kDa Polyacrylonitrile (PAN) Membrane; aqueous; Size: 25 x 282 mm	10

<b>AF2000 MF Series (Asymmetric Flow FFF) Mini Version</b>		<b>Qty / Pack</b>
<b>AF4-Membranes for Aqueous/Organic Solvents</b>		
Z-AF4-MEM-721-300D	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 300 Da cut-off; Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-721-1KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 1 kDa cut-off; Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-721-4KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 4 kDa cut-off; Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-721-5KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 5 kDa cut-off; Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-721-10KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-721-20KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-721-50KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 50 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-721-150KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-722-1KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 1 kDa cut-off; reg. Cellulose (RC), Size 25 x 149 mm	10
Z-AF4-MEM-722-5KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 5 kDa cut-off; reg. Cellulose (RC), Size 25 x 149 mm	10
Z-AF4-MEM-722-10KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 10 kDa cut-off; reg. Cellulose (RC), Size 25 x 149 mm	10
Z-AF4-MEM-722-30KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 30 kDa cut-off; reg. Cellulose (RC), Size 25 x 149 mm	10
Z-AF4-MEM-722-100KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 100 kDa cut-off; reg. Cellulose (RC), Size 25 x 149 mm	10
Z-AF4-MEM-725-10KD	NovaRC AMPH amphiphilic Membrane pre-cut for 700 Channel Series, aqueous, 10kDa, reg.Cellulose, for amphiphilic polymers/cationic polyelectrolytes, Size 25 x 149 mm	10
Z-AF4-MEM-723-10KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Cellulose Triacetate (CTA), Size 25 x 149 mm	10
Z-AF4-MEM-723-20KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off, Cellulose Triacetate (CTA), Size 25 x 149 mm	10
Z-AF4-MEM-724-150KD	NovaTF AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off, PVDF hydrophilic, Size 25 x 149 mm	10
Z-AF4-MEM-727-1KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, 1 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 25 x 149 mm	10
Z-AF4-MEM-727-3KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, 3 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 25 x 149 mm	10
Z-AF4-MEM-727-5KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, 5 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 25 x 149 mm	10
Z-AF4-MEM-727-10KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, 10 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 25 x 149 mm	10
Z-AF4-MEM-727-30KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, 30 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 25 x 149 mm	10

<b>AF2000 MF Series (Asymmetric Flow FFF) Micro Version</b>		<b>Qty / Pack</b>
<b>AF4-Membranes for Aqueous/Organic Solvents</b>		
Z-AF4-MEM-731-10KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Poly(ether)sulfone (PES), Size 75 x 15 mm	10
Z-AF4-MEM-732-5KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 5 kDa cut-off; reg. Cellulose (RC), Size 75 x 15 mm	10
Z-AF4-MEM-732-10KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 10 kDa cut-off; reg. Cellulose (RC), Size 75 x 15 mm	10
Z-AF4-MEM-735-10KD	NovaRC AMPH amphiphilic AQU Membrane pre-cut for 700 Channel Series, aqueous, 10kDa, reg.Cellulose, for amphiphilic polymers/cationic polyelectrolytes, Size 75 x 15 mm	10
Z-AF4-MEM-737-1KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, 1 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 75 x 15 mm	10
Z-AF4-MEM-737-3KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, 3 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 75 x 15 mm	10
Z-AF4-MEM-737-5KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, 5 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 75 x 15 mm	10
Z-AF4-MEM-737-10KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, 10 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 75 x 15 mm	10
Z-AF4-MEM-737-30KD	NovaPAN AQU Membrane pre-cut for 700 Channel Series, 30 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 75 x 15 mm	10