

Membranes for Flow Field-Flow Fractionation

Postnova offers a broad range of special membranes for Flow FFF Systems. Please check the system you use and select the corresponding membrane from the list below or contact our customer support.

AF2000 MF Series (Asymmetric Flow FFF) Micro Version		Qty / Pack
AF4-Membranes for Aqueous/Organic Solvents		
Z-AF4-MEM-731-10KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Poly(ether)sulfone (PES), Size 75 x 15 mm	10
Z-AF4-MEM-732-5KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 5 kDa cut-off; reg. Cellulose (RC), Size 75 x 15 mm	10
Z-AF4-MEM-732-10KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 10 kDa cut-off; reg. Cellulose (RC), Size 75 x 15 mm	10
Z-AF4-MEM-735-10KD	NovaRC AMPH amphiphilic AQU Membrane pre-cut for 700 Channel Series, aqueous, 10kDa, reg.Cellulose, for amphiphilic polymers/cationic polyelectrolytes, Size 75 x 15 mm	10
Z-AF4-MEM-737-1KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, 1 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 75 x 15 mm	10
Z-AF4-MEM-737-3KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, 3 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 75 x 15 mm	10
Z-AF4-MEM-737-5KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, 5 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 75 x 15 mm	10
Z-AF4-MEM-737-10KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, 10 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 75 x 15 mm	10
Z-AF4-MEM-737-30KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, 30 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 75 x 15 mm	10

AF2000 MF Series (Asymmetric Flow FFF) Mini Version		Qty / Pack
AF4-Membranes for Aqueous/Organic Solvents		
Z-AF4-MEM-721-300D	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 300 Da cut-off; Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-721-1KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 1 kDa cut-off; Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-721-4KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 4 kDa cut-off; Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-721-5KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 5 kDa cut-off; Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-721-10KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-721-20KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-721-50KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 50 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-721-150KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 149 mm	10
Z-AF4-MEM-722-1KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 1 kDa cut-off; reg. Cellulose (RC), Size 25 x 149 mm	10
Z-AF4-MEM-722-5KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 5 kDa cut-off; reg. Cellulose (RC), Size 25 x 149 mm	10
Z-AF4-MEM-722-10KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 10 kDa cut-off; reg. Cellulose (RC), Size 25 x 149 mm	10
Z-AF4-MEM-722-30KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 30 kDa cut-off; reg. Cellulose (RC), Size 25 x 149 mm	10
Z-AF4-MEM-722-100KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 100 kDa cut-off; reg. Cellulose (RC), Size 25 x 149 mm	10
Z-AF4-MEM-725-10KD	NovaRC AMPH amphiphilic Membrane pre-cut for 700 Channel Series, aqueous, 10kDa, reg.Cellulose, for amphiphilic polymers/cationic polyelectrolytes, Size 25 x 149 mm	10
Z-AF4-MEM-723-10KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Cellulose Triacetate (CTA), Size 25 x 149 mm	10
Z-AF4-MEM-723-20KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off, Cellulose Triacetate (CTA), Size 25 x 149 mm	10
Z-AF4-MEM-724-150KD	NovaTF AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off, PVDF hydrophilic, Size 25 x 149 mm	10
Z-AF4-MEM-727-1KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, 1 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 25 x 149 mm	10
Z-AF4-MEM-727-3KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, 3 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 25 x 149 mm	10
Z-AF4-MEM-727-5KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, 5 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 25 x 149 mm	10
Z-AF4-MEM-727-10KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, 10 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 25 x 149 mm	10
Z-AF4-MEM-727-30KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, 30 kDa Polyacrylnitrile (PAN) Membrane; Aqueous; Size: 25 x 149 mm	10

AF2000 MF Series (Asymmetric Flow FFF) Analytical Version		Qty / Pack
AF4-Membranes for Aqueous/Organic Solvents		
Z-AF4-MEM-711-300D	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 300 Da cut-off; Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-711-1KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 1 kDa cut-off; Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-711-4KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 4 kDa cut-off; Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-711-5KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 5 kDa cut-off; Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-711-10KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-711-20KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-711-50KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 50 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-711-150KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off, Poly(ether)sulfone (PES), Size 25 x 282 mm	10
Z-AF4-MEM-712-1KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 1 kDa cut-off; reg. Cellulose (RC), Size 25 x 282 mm	10
Z-AF4-MEM-712-5KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 5 kDa cut-off; reg. Cellulose (RC), Size 25 x 282 mm	10
Z-AF4-MEM-712-10KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 10 kDa cut-off; reg. Cellulose (RC), Size 25 x 282 mm	10
Z-AF4-MEM-712-30KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 30 kDa cut-off; reg. Cellulose (RC), Size 25 x 282 mm	10
Z-AF4-MEM-712-100KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 100 kDa cut-off; reg. Cellulose (RC), Size 25 x 282 mm	10
Z-AF4-MEM-715-10KD	NovaRC AMPH amphiphilic Membrane pre-cut for 700 Channel Series, aqueous, 10kDa, reg.Cellulose, for amphiphilic polymers/cationic polyelectrolytes, Size 25 x 282 mm	10
Z-AF4-MEM-713-10KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Cellulose Triacetate (CTA), Size 25 x 282 mm	10
Z-AF4-MEM-713-20KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off, Cellulose Triacetate (CTA), Size 25 x 282 mm	10
Z-AF4-MEM-714-50KD	NovaTF AQU Membrane pre-cut for 700 Channel Series, aqueous, 50 kDa cut-off, PVDF hydrophilic, Size 25 x 282 mm	10
Z-AF4-MEM-714-150KD	NovaTF AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off, PVDF hydrophilic, Size 25 x 282 mm	10
Z-AF4-MEM-717-1KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, 1 kDa Polyacrylonitrile (PAN) Membrane; aqueous; Size: 25 x 282 mm	10
Z-AF4-MEM-717-3KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, 3 kDa Polyacrylonitrile (PAN) Membrane; aqueous; Size: 25 x 282 mm	10
Z-AF4-MEM-717-5KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, 5 kDa Polyacrylonitrile (PAN) Membrane; aqueous; Size: 25 x 282 mm	10
Z-AF4-MEM-717-10KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, 10 kDa Polyacrylonitrile (PAN) Membrane; aqueous; Size: 25 x 282 mm	10
Z-AF4-MEM-717-30KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, 30 kDa Polyacrylonitrile (PAN) Membrane; aqueous; Size: 25 x 282 mm	10

AF2000 MF Series (Asymmetric Flow FFF) Semipreparative Version		Qty / Pack
AF4-Membranes for Aqueous/Organic Solvents		
Z-AF4-MEM-741-300D	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 300 Da cut-off; Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-741-1KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 1 kDa cut-off; Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-741-4KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 4 kDa cut-off; Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-741-5KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 5 kDa cut-off; Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-741-10KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-741-20KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off, Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-741-50KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 50 kDa cut-off, Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-741-150KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off, Poly(ether)sulfone (PES), Size 56 x 283 mm	1
Z-AF4-MEM-742-1KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 1 kDa cut-off; reg. Cellulose (RC), Size 56 x 283 mm	1
Z-AF4-MEM-742-5KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 5 kDa cut-off; reg. Cellulose (RC), Size 56 x 283 mm	1
Z-AF4-MEM-742-10KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 10 kDa cut-off; reg. Cellulose (RC), Size 56 x 283 mm	1
Z-AF4-MEM-742-30KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 30 kDa cut-off; reg. Cellulose (RC), Size 56 x 283 mm	1
Z-AF4-MEM-742-100KD	NovaRC ORG Membrane pre-cut for 700 Channel Series, organic, 100 kDa cut-off; reg. Cellulose (RC), Size 56 x 283 mm	1
Z-AF4-MEM-745-10KD	NovaRC AMPH amphiphilic Membrane pre-cut for 700 Channel Series, aqueous, 10kDa cut-off, reg.Cellulose, for amphiphilic polymers/cationic polyelectrolytes, Size 56 x 283 mm	1
Z-AF4-MEM-743-10KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Cellulose Triacetate (CTA), Size 56 x 283 mm	1
Z-AF4-MEM-743-20KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off, Cellulose Triacetate (CTA), Size 56 x 283 mm	1
Z-AF4-MEM-744-150KD	NovaTF AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off, PVDF hydrophilic, Size 56 x 283 mm	1
Z-AF4-MEM-747-1KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, aqueous, 1 kDa cut-off, Polyacrylnitrile (PAN), Size 56 x 283 mm	1
Z-AF4-MEM-747-3KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, aqueous, 3 kDa cut-off, Polyacrylnitrile (PAN), Size 56 x 283 mm	1
Z-AF4-MEM-747-5KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, aqueous, 5 Da cut-off, Polyacrylnitrile (PAN), Size 56 x 283 mm	1
Z-AF4-MEM-747-10KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Polyacrylnitrile (PAN), Size 56 x 283 mm	1
Z-AF4-MEM-747-30KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, aqueous, 30 kDa cut-off, Polyacrylnitrile (PAN), Size 56 x 283 mm	1

AF2000 MF Series (Asymmetric Flow FFF) Preparative Version		Qty / Pack
AF4-Membranes for Aqueous/Organic Solvents		
Z-AF4-MEM-601-300D	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 300 Da cut-off; Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-601-1KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 1 kDa cut-off; Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-601-4KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 4 kDa cut-off; Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-601-5KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 5 kDa cut-off; Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-601-10KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-601-20KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off, Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-601-50KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 50 kDa cut-off, Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-601-150KD	NovaPES AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off, Poly(ether)sulfone (PES), Size 110 x 294 mm	1
Z-AF4-MEM-602-1KD	NovaRC AQU Membrane pre-cut for 700 Channel Series, aqueous, 1 kDa cut-off; reg. Cellulose (RC), Size 110 x 294 mm	1
Z-AF4-MEM-602-5KD	NovaRC AQU/ORG Membrane pre-cut for 700 Channel Series, aqueous & organic, 5 kDa cut-off; reg. Cellulose (RC), Size 110 x 294 mm	1
Z-AF4-MEM-602-10KD	NovaRC AQU/ORG Membrane pre-cut for 700 Channel Series, aqueous & organic, 10 kDa cut-off; reg. Cellulose (RC), Size 110 x 294 mm	1
Z-AF4-MEM-602-30KD	NovaRC AQU Membrane pre-cut for 700 Channel Series, aqueous, 30 kDa cut-off; reg. Cellulose (RC), Size 110 x 294 mm	1
Z-AF4-MEM-602-100KD	NovaRC AQU Membrane pre-cut for 700 Channel Series, aqueous, 100 kDa cut-off; reg. Cellulose (RC), Size 110 x 294 mm	1
Z-AF4-MEM-605-10KD	NovaRC AMPH amphiphilic Membrane pre-cut for 700 Channel Series, aqueous, 10kDa cut-off, reg.Cellulose, for amphiphilic polymers/cationic polyelectrolytes, Size 110 x 294 mm	1
Z-AF4-MEM-603-10KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Cellulose Triacetate (CTA), Size 110 x 294 mm	1
Z-AF4-MEM-603-20KD	NovaCTA AQU Membrane pre-cut for 700 Channel Series, aqueous, 20 kDa cut-off, Cellulose Triacetate (CTA), Size 110 x 294 mm	1
Z-AF4-MEM-604-150KD	NovaTF AQU Membrane pre-cut for 700 Channel Series, aqueous, 150 kDa cut-off, PVDF hydrophilic, Size 110 x 294 mm	1
Z-AF4-MEM-607-1KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, aqueous, 1 kDa cut-off, Polyacrylnitrile (PAN), Size 110 x 294 mm	1
Z-AF4-MEM-607-3KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, aqueous, 3 kDa cut-off, Polyacrylnitrile (PAN), Size 110 x 294 mm	1
Z-AF4-MEM-607-5KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, aqueous, 5 Da cut-off, Polyacrylnitrile (PAN), Size 110 x 294 mm	1
Z-AF4-MEM-607-10KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, aqueous, 10 kDa cut-off, Polyacrylnitrile (PAN), Size 110 x 29 mm	1
Z-AF4-MEM-607-30KD	NovaMEM AQU Membrane pre-cut for 700 Channel Series, aqueous, 30 kDa cut-off, Polyacrylnitrile (PAN), Size 110 x 294 mm	1